Autonorm™ Human Liquid L-1

REF 202615

LOT 1510483

2017-12

10 x 8 mL

IVD

EN

Product description

Autonorm™ Human Liquid is a liquid human-based control serum produced from blood collected from voluntary blood donors. Each unit is separately checked using CE approved tests to be negative for HBs antigen, HIV p24-antigen and HIV I, II and HCV antibodies. However, since no method can completely exclude the presence of infectious agents, this material should be handled as an ordinary patient sample.

Safety Data Sheet (SDS) is available at our website www.sero.no.

Intended use

Autonorm™ Human Liquid is intended for use as an assayed quality control material to monitor the precision of laboratory measurement procedures.

Instructions for use

To secure uniformity, ensure that each vial is similarly treated with respect to thawing. Always mix contents thoroughly before use.

1. Thaw the vial contents by gentle continuous rolling for approximately 40 minutes until the material has thawed completely and achieved room temperature.
2. Remove the amount of material required for testing.
3. Promptly replace the stopper in the vial and store the remaining control material at 2-8 °C.
4. Any material removed from the vial but not used should be discarded. DO NOT RETURN control material to the original vial once removed.
5. Establish target values and analytical ranges for the control material based on the measurement procedures in the lab.

Storage and stability

This product will be stable until the expiry date when stored unopened at ≤ -20 °C with the following exception: The activity of ALAT, Bilirubin total, CK and LDH may decrease about 1% per month.

Once thawed and opened, all analytes are stable for 14 days when stored tightly capped in the original vial at 2-8 °C. CK is sensitive to temperature changes. CK and Bilirubin are sensitive to light.

Vials aliquoted and refrozen within 30 minutes at ≤ -20 °C in screw cap tubes with O-ring are stable for 1 month. Discard the remaining material after each use.

Limitations

All stability data require that bacterial contamination is avoided. Increased turbidity may indicate bacterial growth.

To avoid Zinc contamination from vial and rubber stopper, transfer the serum to plastic tubes after reconstitution.

Component levels

No precise analytical values have been elaborated for Autonorm™ Human Liquid. It is therefore recommended that each laboratory establishes its own analytical mean and acceptable range and only use the provided values as a guide.

Note!

Variability of HDL and LDL has been observed with some reagent suppliers.

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