

11. Assay Procedure Continued...

- After 20 minutes, decant or aspirate the well contents and wash the wells 3 times using an automatic plate washer or the manual wash procedure (see below). Careful washing is the key to good results. Blot the wells on absorbent paper before proceeding. **Do not allow the wells to dry out.**

Manual Wash Procedure

Empty the wells by inversion. Using a multi-channel pipette or wash bottle, fill the wells with wash buffer. Empty by inversion and blot the wells on absorbent paper. Repeat this wash process 2 more times.

- Dispense 100µl of Conjugate (**Reagent 3**) into each well. Incubate the wells for **20 minutes** at room temperature.
- After 20 minutes, discard the well contents and carefully wash the wells 4 times with Wash Buffer. Ensure that the wells are empty but do not allow to dry out.
- Using a repeating dispenser, rapidly dispense 100µl of TMB Substrate (**Reagent 4**) into each well. Incubate the plate for **10 minutes**.
- Add 100µl of Stop Solution (**Reagent 5**) to each well. To allow equal reaction times, the Stop Solution should be added to the wells in the same order as the TMB Substrate.
- Read the optical density (OD) of each well at 450nm in a microplate reader within 10 minutes. A 620nm filter may be used as a reference wavelength.

12. Quality Control

Quality control data is supplied on the lot-specific QC certificate included in the kit. Controls are intended to monitor for substantial reagent failure. Any well positive by spectrophotometer but without visible colour should be cleaned on the underside and re-read. If OD-values below zero are observed, the wavelengths used should be verified, the reader re-blanked to air and the measurements repeated.

13. Interpretation of Results

Negative samples: OD < OD of 10 U/ml standard
Positive samples: OD >= OD of 10 U/ml standard

- A negative result indicates no current or previous infection with CMV. Such individuals are presumed to be susceptible to primary infection. However, specimens taken too early during a primary infection may not have detectable levels of IgM antibody. If a primary infection is suspected, another specimen should be taken within 7 days and tested concurrently in the same assay with the original specimen to look for seroconversion.
- A positive result indicates a primary or a reactivated infection with CMV. Such individuals are presumed to be at risk of transmitting CMV infection.

14. Limitations of the Procedure

- Results of the Genesis Diagnostics CMV IgM assay are not by themselves diagnostic and should be interpreted in conjunction with the patient's clinical condition and the results of other diagnostic tests.
- A negative result does not rule out a primary or reactivated infection with CMV.
- CMV-specific IgM antibody usually does not develop until the patient has been clinically ill for a week or more. Therefore, samples taken too early in the course of a primary infection may not have detectable levels of IgM.

- Patients may continue to produce CMV specific IgM antibody for 6-9 months following a primary infection.
- Heterotypic IgM antibody responses may occur in patients with Epstein-Barr virus and give false positive results.
- The incidence of CMV infection varies with age, geographical location, sexual behaviour and socio-economic class. CMV specific IgM usually develops after a patient has been clinically ill for at least a week or more. Most patients produce IgM transiently within 16 weeks of seroconversion. However, some patients continue to produce IgM for 6-9 months after seroconversion.

15. Performance Characteristics

Comparative Study

The Genesis Diagnostics CMV IgM kit was compared with another commercially available ELISA procedure for the detection of IgM antibodies to CMV. The Genesis kit showed 100% agreement with the other ELISA. The results are summarised below.

n=70	Reference +	Reference -
Genesis +	10	0
Genesis -	0	60

16. Reproducibility

Within Assay Precision

CV%: <12%

Between Assay Precision

CV%: <12%

17. Method Summary

- Mix IgG absorbent and Sample Diluent 1:4 and dilute all samples 1:100
- Dispense 100µl of the 10 U/ml standard, each control and diluted sample into the microplate wells
- Incubate for **20 minutes** at room temperature.
- Wash the wells three times
- Dispense 100µl of Conjugate (**Reagent 3**) into each well
- Incubate at room temperature for **20 minutes**
- Wash the wells four times
- Add 100µl of TMB Substrate (**Reagent 4**) to each well
- Incubate at room temperature for **10 minutes**
- Add 100µl Stop Solution (**Reagent 5**) to each well
- Read the optical density at 450nm

18. Further Reading

Drew WL: Diagnosis of cytomegalovirus infection. Rev Infect Dis 10:5468-5475, 1988
Booth JC *et al*: Comparison of enzyme-linked immunosorbent assay, radioimmunoassay, complement fixation, anticomplement immunofluorescence and passive agglutination techniques for detecting cytomegalovirus IgG antibody. J Clin Pathol 35: 1345-1348, 1982
Dylewski JS *et al*: Large scale serological screening for cytomegalovirus antibodies in homosexual males by enzyme-linked immunosorbent assay. J Clin Micro 19: 200-203, 1984
Stagno S *et al*: Immuno-globulin M antibodies detected by enzyme-linked immunosorbent assay in the diagnosis of cytomegalovirus infections in pregnant women and new born infants. J Clin Micro 21: 930-935, 1985

Eden Research Park, Henry Crabb Road, Littleport, Cambridgeshire,
CB6 1SE, UK Tel+ 44(0)1353 862220 Fax+44(0)1353 863330
Email: support@elisa.co.uk Web: www.omegadiagnostics.com
Certified to ISO9001:2008, ISO13485:2003

Genesis Diagnostics Ltd is a subsidiary of Omega Diagnostics Group PLC

G • E • N • E • S • I • S
Diagnostics

CMV IgM ELISA Kit

Qualitative assay for anti-CMV IgM antibodies

Product Code: GD085

For Research Use Only. Not for use in diagnostic procedures

Not for sale or use in the EU

1. Materials Included in the Kit

- Microplate** 96 wells in 12 x 8 break-apart strips, pre-coated with density gradient-purified CMV antigen.
- Reagent 1: Sample Diluent** 10mM Tris-buffered saline, pH 7.2 with antimicrobial agent, 46ml, (blue). Read the instructions before use.
- IgG absorbent:** Anti-human IgG, 3 x 3.5ml. Read the instructions before use.
- Reagent 2: Wash Buffer** 100mM Tris-buffered saline with detergent, pH 7.2, 100 ml, **concentrate** (x10)
- Reagent 3: Conjugate** goat anti-human IgM (green) conjugated to horseradish peroxidase in protein stabilising solution and antimicrobial agent, 12ml, ready to use
- Reagent 4: TMB Substrate** aqueous solution of TMB and hydrogen peroxide, 12 ml, ready to use
- Reagent 5: Stop Solution** 0.25M sulphuric acid, 12 ml, ready to use
- Standards:** 10 U/ml (yellow), 1ml of 10mM Tris-buffered saline with human serum IgM antibodies Rubella, ready to use
- Positive control:** 1ml of 10mM Tris-buffered saline containing human serum antibodies to Rubella, 100 U/ml, (red), ready to use.
- Negative control:** 1ml of 10mM Tris-buffered saline containing normal human serum, (green), ready to use.
- Instructions for use**

2. Other Equipment Required

Test tubes for dilution • graduated cylinder for preparing wash buffer • precision pipettes and disposable tips to deliver 5µl, 100µl, 1ml • EIA microplate washer or multi-channel pipette or wash bottle • distilled or de-ionised water • absorbent paper • EIA microplate reader with 450nm and optional 620nm reference filter. Alternatively, a suitable, self-validated automated system may be used.

Instrumentation, whether manual or automated, should meet the following criteria: pipettes with better than 3% imprecision with no carry over between pipetting steps; microplate washers should remove 99% of fluid; automated machines should minimise time between washing and adding the next reagent.

3. Intended Use

The CMV IgM kit is a rapid ELISA designed for the qualitative detection of IgM antibodies to cytomegalovirus (CMV) in human serum or plasma. The assay is intended to be used to evaluate serologic evidence of primary or reactivated infection with CMV, and is for *in vitro* research use only.

4. Explanation of the Test

CMV infections can be classified as congenital (acquired before birth), perinatal (acquired at birth) or postnatal (acquired after birth). 95% of newborn infants congenitally infected with CMV exhibit no clinically overt disease at birth. Disease in the remaining 5% can be severe and may result in neurological damage in survivors or death early after birth.

The prognosis for congenitally infected infants who are asymptomatic at birth is variable. Many subsequently develop hearing loss and varying degrees of mental retardation and central nervous system disorders. Surveys show the incidence congenital CMV infection to be between 0.5 – 2.5%.

Perinatally infected infants start shedding the virus 3–12 weeks after delivery and generally remain asymptomatic. Postnatal CMV infections are usually asymptomatic and are acquired through close contact with individuals who are shedding the virus. A small percentage of individuals develop a negative heterophile-antibody infectious mononucleosis syndrome characterised by fever, lethargy and atypical lymphocytosis.

In immunocompromised patients e.g. allograft recipients, cancer patients and AIDS patients, CMV infections occur frequently, often from re-activation of latent infection, and may be life-threatening. CMV infections may occur following blood transfusions. Most transfusion acquired CMV infections are either sub-clinical or characterised by CMV mononucleosis. However, in immunocompromised patients, considerable morbidity and mortality can result from a transfusion-acquired CMV infection.

Serologic tests which measure IgG antibodies to CMV can aid in the diagnosis of CMV infection when paired acute and convalescent sera are tested simultaneously and seroconversion or a significant rise in titre can be demonstrated. Also serologic testing of blood transfusion donors and recipients can help prevent transfusion-acquired infection.

IgM antibodies are produced during the first 2 to 3 weeks of infection with CMV and exist only transiently in most individuals. Serologic procedures that measure the presence of IgM antibodies help to discriminate between primary

and recurrent infections since IgM antibodies are rarely found in recurrent infection.

5. Principle of the Test

Test sera are diluted (1:100) with the sample diluent provided. Anti-human IgG is added to the sample diluent sample to eliminate the possibility of interference by antigen-specific IgG and rheumatoid factor, if present. Diluted serum or plasma specimens are incubated for 20 minutes to allow specific antibodies to CMV to bind to the antigen-coated wells. After washing away unbound antibodies and other serum constituents, CMV specific IgM is detected using rabbit anti-human IgM conjugated to horseradish peroxidase. After 20 minutes incubation, unbound conjugate is removed by washing, and TMB enzyme substrate is added for 10 minutes. A blue colour develops if antibodies to CMV are present. Addition of stop solution gives a yellow colour and the optical densities of controls, 10 U/ml standard and samples are measured using a microplate reader.

6. Safety Precautions

1. Only experienced laboratory personnel should use this test. The test protocol must be followed strictly.
2. CAUTION: the device contains material of human and animal origin and should be handled as a potential transmitter of diseases. All human source material used in the preparation of standards and control for this product have been tested and found negative by ELISA for antibodies to HIV, HbsAg and HCV. No test method, however, can offer complete assurance that infectious agents are absent. Therefore, all reagents containing human material should be handled as if potentially infectious. Operators should wear gloves and protective clothing when handling any patient sera or serum based products.
3. Reagents of this kit contain antimicrobial agents and the Substrate solution contains 3,3',5,5'-tetramethylbenzidine. Avoid contact with the skin and eyes. Rinse immediately with plenty of water if any contact occurs.
4. The Stop Solution contains 0.25M sulphuric acid. Avoid contact with skin and eyes. Rinse immediately with plenty of water if contact occurs.
5. Any liquid that has been brought into contact with potentially infectious material has to be discarded in a container with a disinfectant. Dispose of plates and specimens as clinical waste. Any unused reagents should be flushed away with copious amounts of water. Disposal must be performed in accordance with local legislation.

7. Technical Precautions

1. Strips and solutions should not be used if the foil bag is damaged or liquids have leaked.
2. Allow all reagents and the microplate to reach room temperature before use. Ensure that the microplate foil bag containing any unused strips is well sealed and contains the desiccant to avoid moisture. Store at 2 – 8°C after use.
3. Include the Positive and Negative Control in every test run to monitor for reagent stability and correct assay performance.
4. Strictly observe the indicated incubation times and temperature.
5. When automating, consider excess volumes required for setting up the instrument and dead volume of robot pipette
6. Ensure that no cross-contamination occurs between wells. Keep all pipettes and other equipment used for Conjugate completely separate from the TMB Substrate reagent.
7. When pipetting Conjugate or TMB Substrate, aliquots for the required numbers of wells should be taken to avoid multiple entry of pipette tips

into the reagent bottles. Never pour unused reagents back into the original bottles.

8. Do not allow microwells to dry between incubation steps.
9. Strictly follow the described wash procedure. Insufficient washing may cause high background signal.
10. Avoid direct sunlight and exposure to heat sources during all incubation steps.
11. Replace colour-coded caps on their correct vials to avoid cross-contamination
12. It is important to dispense all samples and controls into the wells without delay. Therefore ensure that all samples are ready to dispense.

8. Shelf Life and Storage Conditions

On arrival, store the kit at 2 - 8°C. Once opened the kit is stable for 3 months (or until its expiry date if less than 3 months). Do not use kits beyond their expiry date. Do not freeze any kit component. The diluted Wash Buffer has a shelf life of 3 months if stored in a closed bottle at 2 - 8°C.

9. Specimen Collection and Storage

Serum and plasma samples may be used and should be stored at -20°C for long-term storage. Frozen samples must be mixed well after thawing and prior to testing. Repeated freezing and thawing can affect results. Addition of preservatives to the serum sample may adversely affect the results. Microbially contaminated, heat-treated or specimens containing particulate matter should not be used. Grossly haemolysed, icteric or lipaemic specimens should be avoided.

10. Preparation of Reagents

1. Prepare only sufficient IgG-absorbent-containing sample diluent for the number of samples to be tested. Add one part IgG absorbent to 4 parts of Sample Diluent (**Reagent 1**) as shown in the examples below and mix thoroughly. Discard any unused IgG-absorbent-containing diluent.

Approx. # of Samples	Vol. of Sample Diluent (ml)	Vol. of IgG Absorbent (ml)
24	10	2.5
48	20	5.0
72	30	7.5
96	40	10.0

2. Dilute the Wash Buffer (**Reagent 2**) 1:9 in distilled water to make sufficient buffer for the assay run e.g. add 50ml wash buffer concentrate to 450ml water.

11. Assay Procedure

1. Dilute patient samples 1:100 (e.g. 5µl serum plus 0.5 ml IgG-absorbent-containing sample diluent). It is important to dispense all samples, standards and controls into the wells without delay. Therefore ensure that all samples are ready to dispense.
2. Assemble the number of strips required for the assay.
3. For qualitative assays, dispense 100 µl of the negative control, the 10 U/ml standard, the positive control and the diluted patient sample into the wells.
4. Incubate for **20 minutes** at room temperature. During all incubations, avoid direct sunlight and close proximity to any heat sources.