N.B. The calibrator (Standard #) values are found in the quality control certificate (SIN283.QC). Calibrator values on **Page 1** should be entered in ascending order, i.e. the lowest value first. **IMPORTANT:** the analyser will only update the calibrator values providing the **Auto Fill** button is pressed after typing the value for calibrator 5 on Page: 1. The assay utilises an extrapolated calibration curve therefore the Auto Rerun Range (Result) (Page 3) and Out-of-Range (Page 4) tables must be manually updated using values supplied in the Quality Control certificate (SIN283.QC)

8.4.2 Control volumes

Transfer 150µL of control fluid into a sample cup and place onto the rack.

8.5 Measuring range

The range of this assay is approximately 30 - 1007 IU/mL. In the event of a rerun the upper limit of the assay range is extended to 1500 IU/mL. These values are dependent on the lot specific value of the calibrator in use.

SPAPLUS analyser dilution	Approximate measuring range (IU/mL)
1/1	30 – 1007
1/10	300 - 1500

9 QUALITY CONTROL

The controls provided should be included in all assay runs. The IgE concentration is stated on the accompanying Quality Control certificate (SIN283.QC). Sample results obtained should only be accepted if the control results are within $\pm 20\%$ of the concentration(s) stated.

10 LIMITATIONS

- 10.1 These kits are not suitable for the measurement of lipaemic or haemolysed samples due to the unpredictable degree of non-specific scatter these sample types may generate. Unexpected results should be confirmed using an alternative assay method.
- 10.2 Customers are strongly advised to run all controls with every batch of samples being assayed. Should a control value be out of range against a stored curve, it is recommended that the control is re-assayed using the same calibration curve. If the control value is still out of range the curve should be recalibrated and the controls re-assayed. If control values are out of range against the new calibration curve check the instrument and parameters entered before repeating the assay. If problems persist, refer to the supplier.
- 10.3 Diagnosis cannot be made and treatment must not be initiated on the basis of IgE measurements alone. Clinical history and other laboratory findings must also be taken into account.
- 10.4 Potential occurrences of antigen excess cannot be completely excluded; in rare cases samples with monoclonal IgE present may give falsely low results due to antigen excess. Where this is possible or suspected, it is recommended that the sample is re-assayed at a higher dilution to confirm the result.

11 EXPECTED VALUES

REFERENCE RANGE^{3,4}

The concentration of IgE in non-atopical subject's serum varies according to age.

Age group	Upper limit of the normal range (95th percentile)
New born	1.5 IU/mL
Up to 1 year	15.0 IU/mL
Children 1 - 5 years	60.0 IU/mL
Children 6 - 9 years	90.0 IU/mL
Children 10 - 15 years	200.0 IU/mL
Adults	100 IU/mL

It is recommended that each laboratory establish its own reference range to reflect the age, sex, diet and geographical location of the population.

For diagnostic purposes IgE results should be assessed in conjunction with the patient's medical history, clinical examination and other findings.

12 PERFORMANCE CHARACTERISTICS

12.1 Precision

IgE precision summary						
	n	Mean (IU/mL)	Intra-assay		Total precision	
			SD	CV %	SD	CV %
Serum 1	80	115.08	3.83	3.3%	6.25	5.4%
Serum 2	80	378.64	8.00	2.1%	11.56	3.1%
Serum 3	80	497.14	9.73	2.0%	14.56	2.9%

12.2 Comparison

This method (Y) was compared with another commercially available method (X) and the following Passing & Bablok equation obtained:

$$y = 1.01x - 8.49$$

and a correlation coefficient of $r = 0.998$

73 patient samples were analysed spanning the range 31.1 to 932.8 IU/mL.

12.3 Analytical sensitivity

The Limit of Quantitation (LoQ), the Limit of Detection (LoD) and the Limit of Blank (LoB) were determined consistent with CLSI guidelines EP17-A2. LoQ is the smallest concentration that can be detected reliably. LoD is the smallest concentration that can be detected to determine the presence or absence of IgE. LoB is the highest concentration that is likely to be observed in a blank sample.

	IgE concentration IU/mL
Limit of Blank	4.92
Limit of Detection	10.29
Limit of Quantitation	30

12.4 Linearity

The method is linear up to an IgE concentration of 911.3 IU/mL.

This gave a regression plot of:

$$y = 1.01x + 3.17; r = 0.995$$

(y = measured IgE concentration, x = theoretical concentration)

12.5 Interference

The analytes below were tested up to the following levels and were found not to interfere:

Haemoglobin 1000 mg/dL
 Intralipid 1000 mg/dL
 Free Bilirubin 60 mg/dL
 Conjugated Bilirubin 60 mg/dL
 Triglycerides 2000 mg/dL

12.6 Antigen excess

Antigen excess effects are not noted until levels approach 16000 IU/mL.

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Nur zur *in-vitro* Diagnostic

Bestell-Nr.: LK014.S

The Binding Site Group Ltd, 8 Calthorpe Road, Edgbaston, Birmingham, B15 1QT, U.K.
www.bindingsite.co.uk

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Item Name 47 IgE	
DATA INFORMATION	
Units	IU/mL
Decimals	2
ANALYSIS	
Type	End ▼
Main W.Length 1	570 ▼
Sub W.Length	800 ▼
CORR.	
Y =	1 X + 0
CALIBRATION	
Type	Logit 3 ▼
Standard	Auto Fill
1 #	4 #
2 #	5 #
3 #	6 #
NORMAL RANGE	
Method	MALE FEMALE
	LOW HIGH LOW HIGH
Serum	[] [] [] []
Urine	[] [] [] []
Plasma	[] [] [] []
CSF	[] [] [] []
Dialysis	[] [] [] []
Other	[] [] [] []
Page : 1	Print Hard Copy Next Page Save Return

Item Name 47 IgE	
DATA PROCESS	
ASPIRATION	
KIND	Single Double
VOLUME	5
SAMPLE REAGENT1 VOL	200 µL
REAGENT2 VOL	100
FACTOR	
Blank correction	Reaction Check
ENDPOINT LIMIT	2
LINEAR CHECK (%)	0
Third mix	OFF ON
Blank	Water - Blank
DILUTION	
Diluent	99: Dil 1 100: Dil 2
Pre Dilution Rate	10
Auto Rerun Dilution Rate High	10
Auto Rerun Dilution Rate Low	10
MONITOR	
0 LEVEL SPAN 1	3
SPAN	3
PROZONE CHECK	
FIRST	[] [] [] []
SECOND	[] [] [] []
THIRD	[] [] [] []
Min dOD [0]	
Page : 2 Print Prev Page Next Page Save Return	

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Item Name 47 IgE	
Auto Rerun SW	
On Off	
Auto Rerun Range (Result)	
On Off On Off	
Lower Higher	
Serum	# #
Urine	# #
Plasma	# #
CSF	# #
Dialysis	# #
Other	# #
Bottle Size (ml)	
24 Items	36 Items
Reagent1 60	Reagent 1 0
Reagent2 R1 21.0	Reagent2 R1 0
Reagent2 R2 11.0	Reagent2 R2 0
Auto Rerun Condition (Absorbance)	
Absorbance Range	
Lower On Off	
Higher On Off	
Prozone Range	
On Off	
Page : 3 Print Prev Page Next Page Save Return	

Item Name 47 IgE						
Out-of-Range Table						
	NEAT		Pre Dilution		Auto-rerun Dilution (*10)	
	ABOVE	BELOW	ABOVE	BELOW	ABOVE	BELOW
	Cal#1	Cal#2-Cal#6	Cal#1	Cal#2-Cal#6	Cal#1	Cal#2-Cal#6
Serum	#	#	#	#	#	#
Urine	#	#	#	#	#	#
Plasma	#	#	#	#	#	#
CSF	#	#	#	#	#	#
Dialysis	#	#	#	#	#	#
Other	#	#	#	#	#	#
Page : 4	Prev Page		Save		Return	

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