Ferritin Kit for use on SPAPLUS®

For *in-vitro* diagnostic use

Product code: LK055.S

The Binding Site Group Ltd, 8 Calthorpe Road, Edgbaston, Birmingham, B15 1QT, U.K.
www.bindingsite.co.uk
Telephone : +44 (0)121 456 9500
Fax : +44 (0)121 456 9749
E-mail: info@bindingsite.co.uk

In Europe and the USA, SPAPLUS® is a registered trademark of The Binding Site Group Ltd, Birmingham, UK.

1 INTENDED USE

For the quantitative *in vitro* determination of ferritin in serum or plasma. This product is suitable for use on the SPAPLUS analyser.

2 SUMMARY AND EXPLANATION

Ferritin originates in the reticuloendothelial system. It consists of a protein shell, apoferritin (MW 445000) containing varying amounts of iron in its core as ferric hydroxide - phosphate complexes. All complexes contain 2 separate subunits, the acidic H (heavy) type and the weakly basic L (light) type subunit. While the basic isoferritins are concerned with the long-term storage of iron, they are present in the liver, spleen and bone marrow, the acidic isoferritins can be found mainly in the placenta, tumour tissue and myocardium.

Ferritin determination is a good indicator of iron stores in both normal persons and in individuals with iron deficiency or those who are undergoing iron therapy. This assay can be used to differentiate hypochromic anaemia (thalassaemia or chronic infection and tumour anaemias).^{1,4}

3 PRINCIPLE

This method is based on the reaction of a sample containing human ferritin and specific antiserum to form an insoluble complex which can be measured turbidimetrically. A calibration curve is constructed from which the concentration of ferritin is determined.

4 REAGENTS

- 4.1 **R1. Ferritin Reaction Buffer** Containing Glycine buffer 170mmol/L pH8.3, Sodium Chloride 100mmol/L.
- 4.2 **R2. Ferritin Reagent** Latex particles coated with anti-ferritin 0.07% w/v, containing Glycine buffer 170mmol/L pH7.3, Sodium Chloride 100mmol/L.
- 4.3 **Ferritin Controls** Supplied at 2 levels, Low and High. Target values and ranges are supplied in the Quality Control certificate. Supplied ready for use.
- 4.4 **Ferritin Calibrator 1-5** Values have been assigned using an immunoturbidimetric method standardized to the WHO International Standard NIBSC 94/572.

5 CAUTION

For *in vitro* diagnostic use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.

Solutions R1 and R2 contain sodium azide. Avoid ingestion or contact with skin or mucous membranes. In the event of skin contact, flush affected area with copious amounts of water. In the event of contact with eyes or if ingested, seek immediate medical attention.

Sodium azide reacts with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents, flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10% sodium hydroxide.

Health and Safety Data Sheets are available on request.

The reagents must be used only for the purpose intended by suitably qualified laboratory personnel, under appropriate laboratory conditions.

6 STORAGE AND STABILITY

R1. Ferritin Reaction Buffer

Supplied ready for use. Stable up to expiry date when stored at 2-8°C. The *Ferritin Reaction Buffer* may be stored, uncapped, on the analyser for up to 30 days, provided that the main power switch (located at the rear of the left hand panel) is left switched on.

R2. Ferritin Reagent

Supplied ready for use. Stable up to expiry date when stored at 2-8°C protected from light. Invert several times before use, avoiding the formation of foam. The *Ferritin Reagent* may be stored, uncapped, on the analyser for up to 30 days, provided that the main power switch (located at the rear of the left hand panel) is left switched on.

Ferritin Controls

OPENED: Store refrigerated (2-8°C). Ferritin control material is stable for 30 days at 2-8°C if kept capped in original container and free from contamination. Only the required amount of product should be removed. After use, any residual product should NOT BE RETURNED to the original vial.

UNOPENED: Store refrigerated (2-8°C). Stable to expiration date printed on individual vials.

Ferritin Calibrator 1-5

The material is stable to expiry when stored at 2-8°C in the original capped containers. Opened vials are stable for 30 days in the absence of bacterial contamination. Only the required amount of product should be removed. After use, any residual product should NOT BE RETURNED to the original vial.

7 SPECIMEN COLLECTION AND PREPARATION

Serum or EDTA plasma: Stable for 7 days when stored at either 15-25°C or 2-8°C. Sample can be frozen once and stored at -15°C to -25°C for 1 year⁵.

8 METHODOLOGY

8.1 Materials provided

- 8.1.1 1 x 100 tests *Ferritin Reagent SPAPLUS*.
- 8.1.2 1 x 100 tests *Ferritin Reaction Buffer SPAPLUS*.
- 8.1.3 1 x *Ferritin SPAPLUS Calibrator 1-5* (5 x 1.0mL calibrator set for ferritin).
- 8.1.4 1 x 1.0mL *Ferritin SPAPLUS Low Control*
- 8.1.5 1 x 1.0mL *Ferritin SPAPLUS High Control*

8.2 Materials required but not provided

- 8.2.1 Equipment for collection and preparation of test samples e.g. sample tubes, centrifuge etc.
- 8.2.2 A fully operational and equipped SPAPLUS analyser.
- 8.2.3 Current analyser operating instructions: SPAPLUS Reference Guide, Insert Code FIN012.
- 8.2.4 Sample Diluent (99: Dil 1). Binding Site Product Code: SN080.S.

8.3 Reagent preparation

Before loading, gently mix by inversion ensuring no foam or bubbles are generated or remain on the surface as these may interfere with reagent aspiration.

8.4 Test procedure

The user should be familiar with the operation of the SPAPLUS analyser before attempting to carry out the test procedures. The analyser should be prepared for use according to the manufacturer's instructions and the assay protocol entered as described below.

For full details of analyser operation refer to the SPAPLUS Reference Guide (FIN012) supplied with the analyser.

8.4.1 Test parameters

Assay parameters are entered into item number 48.

Item Name 48 FERR		DATA INFORMATION		CALIBRATION	
Units	ng/mL	Type	Logit 2		<input type="checkbox"/> Auto Fill
Decimals	2	Standard			
ANALYSIS		1 #	4 #		
Type	End	2 #	5 #		
		3 #	6 #		
Main W.Length	570	NORMAL RANGE			
Sub W.Length			MALE	FEMALE	
Method		LOW	HIGH	LOW	HIGH
CORR.		Serum	[]	[]	[]
		Urine	[]	[]	[]
		Plasma	[]	[]	[]
		CSF	[]	[]	[]
		Dialysis	[]	[]	[]
		Other	[]	[]	[]
Y =	1 X + 0	Page : 1 Print Hard Copy Next Page Save Return			

Item Name 48 FERR		DATA PROCESS		ABSORBANCE LIMIT	
ASPIRATION		START	END		
KIND	<input type="radio"/> Single <input checked="" type="radio"/> Double	MAIN	51 52	LOW	-3
		SUB	33 34	HIGH	3
SAMPLE	12	FACTOR			
REAGENT1 VOL	140 µL	Blank correction	<input type="radio"/> ON <input checked="" type="radio"/> OFF	Reaction Check	
REAGENT2 VOL	70	ENDPOINT LIMIT	2	LOW	-3
		LINEAR CHECK (%)	0	HIGH	3
Third mix	<input checked="" type="radio"/> OFF <input type="radio"/> ON	DILUTION			
Blank	<input checked="" type="radio"/> Water - Blank	Diluent	<input type="radio"/> 99: Dil 1 <input type="radio"/> 100: Dil 2		
		Pre Dilution Rate			
		Auto Rerun Dilution Rate High	10		
		Auto Rerun Dilution Rate Low			
MONITOR		PROZONE CHECK			
0 LEVEL SPAN	1	START	END	LIMIT (%)	Min dOD [0]
SPAN	3	FIRST	[] [] []		
		SECOND	[] [] []	<input type="radio"/> Low <input checked="" type="radio"/> High	
		THIRD	[] [] []	<input type="radio"/> Low <input checked="" type="radio"/> High	
Page : 2 Print		Prev Page Next Page Save Return			

*Automatically calculated

Item Name 48 FERR		Auto Rerun SW		Auto Rerun Condition (Absorbance)	
		<input checked="" type="radio"/> On <input type="radio"/> Off			
Auto Rerun Range (Result)		Absorbance Range			
		<input type="radio"/> On <input checked="" type="radio"/> Off <input type="radio"/> On <input type="radio"/> Off		Lower	<input type="radio"/> On <input checked="" type="radio"/> Off
		Lower	Higher	Higher	<input type="radio"/> On <input checked="" type="radio"/> Off
Serum	#	#	#		
Urine	#	#	#		
Plasma	#	#	#		
CSF	#	#	#		
Dialysis	#	#	#		
Other	#	#	#		
Bottle Size (ml)		Prozone Range			
24 Items		36 Items		<input type="radio"/> On <input checked="" type="radio"/> Off	
Reagent1	60	Reagent 1	0		
Reagent2 R1	15.0	Reagent2 R1	0		
Reagent2 R2	8.0	Reagent2 R2	0		
Page : 3 Print		Prev Page Save Return			

Item Name 48 FERR						
Out-of-Range Table						
	NEAT		Pre Dilution		Auto-rerun Dilution (*10)	
	ABOVE	BELOW	ABOVE	BELOW	ABOVE	BELOW
	Cal#1	Cal#2-Cal#6	Cal#1	Cal#2-Cal#6	Cal#1	Cal#2-Cal#6
Serum	#	#	#	#	#	#
Urine	#	#	#	#	#	#
Plasma	#	#	#	#	#	#
CSF	#	#	#	#	#	#
Dialysis	#	#	#	#	#	#
Other	#	#	#	#	#	#
Page : 4	Prev Page		Save		Return	

N.B. The calibrator (Standard #) values are found in the Quality Control certificate (SIN284.QC). Calibrator values on Page 1 should be entered in ascending order, i.e. the lowest value first. IMPORTANT: the analyser will only update the calibrator values providing the **Auto Fill** button is pressed after typing the value for calibrator 5 on Page: 1. The assay utilises an extrapolated calibration curve therefore the Auto Rerun Range (Result) (Page 3) and Out-of-Range (Page 4) tables must be manually updated using values supplied in the Quality Control certificate (SIN284.QC)

8.4.2 Control volumes

Transfer 150µL of control fluid into a sample cup and place onto the rack.

8.5 Measuring range

The range of this assay is approximately 7 - 434 ng/mL. In the event of a rerun, the upper limit of the assay range is increased to 868 ng/mL. These values are dependent on the lot specific values of the calibrators in use.

SPAPLUS analyser dilution	Approximate measuring range (ng/mL)
1/1	7 - 434
1/10	70 - 868

9 QUALITY CONTROL

The controls provided should be included in all assay runs. The ferritin concentration is stated on the accompanying Quality Control certificate (SIN284.QC). Sample results obtained should only be accepted if the control results are within ±20% of the concentration(s) stated.

10 LIMITATIONS

- 10.1 These kits are not suitable for the measurement of lipaemic or haemolysed samples due to the unpredictable degree of non-specific scatter these sample types may generate. Unexpected results should be confirmed using an alternative assay method.
- 10.2 Customers are strongly advised to run all controls with every batch of samples being assayed. Should a control value be out of range against a stored curve, it is recommended that the control is re-assayed using the same calibration curve. If the control value is still out of range the curve should be recalibrated and the controls re-assayed. If control values are out of range against the new calibration curve check the instrument and parameters entered before repeating the assay. If problems persist, refer to the supplier.
- 10.3 Diagnosis cannot be made and treatment must not be initiated on the basis of ferritin measurements alone. Clinical history and other laboratory findings must also be taken into account.

11 EXPECTED VALUES

Normal values^{5,6}

Adult Male	20 – 300 ng/mL
Adult Female	10 – 120 ng/mL

A value of less than 12 ng/mL indicates iron deficiency anaemia.

It is recommended that each laboratory establish its own reference range to reflect the age, sex, diet and geographical location of the population.

12 PERFORMANCE CHARACTERISTICS

12.1 Precision

Ferritin precision summary						
	n	Mean (ng/mL)	Intra assay		Total precision	
			SD	CV %	SD	CV %
Serum 1	80	66.67	1.95	2.9%	2.32	3.5%
Serum 2	80	140.97	1.58	1.1%	2.07	1.5%
Serum 3	80	322.74	4.41	1.4%	4.77	1.5%

12.2 Comparison

This method (Y) was compared with another commercially available method (X) and the following Passing & Bablok equation obtained:

$$y = 0.96x + 0.60$$

and a correlation coefficient of $r = 0.996$

95 patient samples were analysed spanning the range 14.1 to 441.6 ng/mL.

12.3 Analytical sensitivity

The limit of Quantitation (LoQ), the limit of Detection (LoD) and the limit of Blank (LoB) were determined consistent with CLSI guidelines EP17-A2. LoQ is the smallest concentration that can be detected reliably. LoD is the smallest concentration that can be detected to determine the presence or absence of ferritin. LoB is the highest concentration that is likely to be observed in a blank sample.

SPAPLUS

	Ferritin concentration ng/mL
Limit of Blank	0.92
Limit of Detection	2.85
Limit of Quantitation	7.00

12.4 Linearity

The method is linear up to a ferritin concentration of 436 ng/mL.

This gave a regression plot of $y = 0.97x + 8.76$; $r = 0.998$ (y = measured ferritin concentration, x = theoretical concentration).

12.5 Interference

The analytes below were tested up to the following levels and were found not to interfere:

Triglyceride 1500 mg/dL
 Haemoglobin 1000 mg/dL
 Free Bilirubin 60 mg/dL
 Conjugated Bilirubin 30 mg/dL
 Intralipid 250 mg/dL

12.6 Antigen excess

Antigen excess effects are not noted until levels approach 30 000 ng/mL.

13 BIBLIOGRAPHY

1. Wick M, Pinggera W, Lehmann, P: Iron metabolism. Diagnosis and Therapy of Anemias, 3rd edition. Vienna/New York : Springer - Verlag, 1996.
2. Albertini A, Arosio P, Chianaone E, Drysdale J (editors). Ferritins and iso-ferritins as biochemical markers. Amsterdam/New York/Oxford Elsevier, 1984.
3. Finlayson NDC. Hereditary (primary) haemochromatosis. BMJ 1990; 301 : 350 - 351.
4. Jacobs D, Kasten B, Demott W, and Wolfson W; Laboratory Test Handbook, 2nd Edition. Williams and Wilkins, 1990.
5. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication. WHO/DIL/LAB/99.1 Rev 2. 2002.
6. Fairbanks, V.F. and Klee, C.G., Textbook of Clinical Chemistry, Tietz, N. W., (ed.), W.B. Saunder Co., Philadelphia, pp 1495-1588 (1986).
7. Cook, J.D., Skikne, B.S., Lynch, S.R and Reusser, M.E., Blood, 68, 726 (1986).

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Zur *in-vitro* Diagnostik

Bestell-Nr.: LK055.S

The Binding Site Group Ltd., 8 Calthorpe Road, Edgbaston, Birmingham, B15 1QT, UK.
www.bindingsite.co.uk

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Item Name 48 FERR		CALIBRATION	
DATA INFORMATION		Type	Logit 2 ▼
Units	ng/mL	Standard	Auto Fill
Decimals	2	1 #	4 #
ANALYSIS		2 #	5 #
Type	End ▼	3 #	6 #
Main W.Length	570 ▼	NORMAL RANGE	
Sub W.Length	▼	MALE FEMALE	
Method		LOW HIGH	LOW HIGH
CORR.		Serum [] []	[] []
		Urine [] []	[] []
		Plasma [] []	[] []
		CSF [] []	[] []
		Dialysis [] []	[] []
		Other [] []	[] []
Y =	1 X + 0	Page : 1 Print Hard Copy Next Page Save Return	

Item Name 48 FERR		DATA PROCESS	
ASPIRATION		READ	ABSORBANCE LIMIT
KIND	○ Single ● Double	MAIN	51 52 LOW -3
		SUB	33 34 HIGH 3
SAMPLE VOLUME		FACTOR	
REAGENT1 VOL	140 µL	Blank correction	* ○ ON ● OFF
REAGENT2 VOL	70 µL	ENDPOINT LIMIT	2 CHECK POINT
		LINEAR CHECK (%)	0 LOW -3
			HIGH 3
Third mix	● OFF ○ ON	DILUTION	
Blank	● Water - Blank	Diluent	● 99: Dil 1 ○ 100: Dil 2
		Pre Dilution Rate	▼
		Auto Rerun Dilution Rate High	10 ▼
		Auto Rerun Dilution Rate Low	▼
MONITOR		PROZONE CHECK	
0 LEVEL SPAN	1	START	END LIMIT (%) Min dOD [0]
SPAN	3	FIRST	[] []
		SECOND	[] [] [] [] ○Low ●High
		THIRD	[] [] [] [] ○Low ●High
Page : 2 Print		Prev Page Next Page Save Return	

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Item Name 48 FERR		Auto Rerun SW	
		Auto Rerun Condition (Absorbance)	
		Absorbance Range	
		Lower	○ On ● Off
		Higher	● On ○ Off
		Prozone Range	
		○ On ● Off	
Bottle Size (ml)			
24 Items	36 Items		
Reagent1	60	Reagent 1	0
Reagent2 R1	15.0	Reagent2 R1	0
Reagent2 R2	8.0	Reagent2 R2	0
Page : 3 Print		Prev Page Save Return	

Item Name 48 FERR						
Out-of-Range Table						
	NEAT		Pre Dilution		Auto-rerun Dilution (*10)	
	ABOVE	BELOW	ABOVE	BELOW	ABOVE	BELOW
	Cal#1	Cal#2-Cal#6	Cal#1	Cal#2-Cal#6	Cal#1	Cal#2-Cal#6
Serum	#	#	#	#	#	#
Urine	#	#	#	#	#	#
Plasma	#	#	#	#	#	#
CSF	#	#	#	#	#	#
Dialysis	#	#	#	#	#	#
Other	#	#	#	#	#	#
Page : 4			Prev Page	Save	Return	

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13 REFERENZEN

- Wick M, Pinggera W, Lehmann, P: Iron metabolism. Diagnosis and Therapy of Anemias, 3rd edition. Vienna/New York : Springer - Verlag, 1996.
- Albertini A, Arosio P, Chianaone E, Drysdale J (editors). Ferritins and iso-ferritins as biochemical markers. Amsterdam/New York/Oxford Elsevier, 1984.
- Finlayson NDC. Hereditary (primary) haemochromatosis. BMJ 1990; 301 : 350 - 351.
- Jacobs D, Kasten B, Demott W, and Wolfson W; Laboratory Test Handbook, 2nd Edition. Williams and Wilkins, 1990.
- Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication. WHO/DIL/LAB/99.1 Rev 2. 2002.
- Fairbanks, V.F. and Klee, C.G., Textbook of Clinical Chemistry, Tietz, N. W., (ed.), W.B. Saunders Co., Philadelphia, pp 1495-1588 (1986).
- Cook, J.D., Skikne, B.S., Lynch, S.R and Reusser, M.E., Blood, 68, 726 (1986).

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The Binding Site Group Ltd, 8 Calthorpe Road, Edgbaston, Birmingham, B15 1QT, U.K.
www.bindingsite.co.uk

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Item Name 48 FERR		CALIBRATION	
DATA INFORMATION		Type Logit 2	Auto Fill
Units	ng/mL	Standard	
Decimals	2	1 #	4 #
ANALYSIS		2 #	5 #
Type	End	3 #	6
Main W.Length 1	570	NORMAL RANGE	
Sub W.Length		MALE	FEMALE
Method		LOW	HIGH
CORR.		Serum	[] [] [] []
Y =		Urine	[] [] [] []
SLOPE		Plasma	[] [] [] []
INTER		CSF	[] [] [] []
1 X + 0		Dialysis	[] [] [] []
		Other	[] [] [] []
Page : 1	Print	Hard Copy	Next Page Save Return

Item Name 48 FERR		DATA PROCESS	
ASPIRATION		READ	ABSORBANCE LIMIT
KIND	Single Double	MAIN	51 52 LOW -3
VOLUME		SUB	33 34 HIGH 3
SAMPLE	12	FACTOR	
REAGENT1 VOL	140	Blank correction	* ON OFF
REAGENT2 VOL	70	ENDPOINT LIMIT	2 CHECK POINT
Third mix		LINEAR CHECK (%)	0 LOW HIGH -3
Blank		DILUTION	
Water - Blank		Diluent	99: Dil 1 100: Dil 2
		Pre Dilution Rate	
		Auto Rerun Dilution Rate High	10
		Auto Rerun Dilution Rate Low	
MONITOR		PROZONE CHECK	
0 LEVEL SPAN 1		FIRST	[] [] [] []
SPAN 3		SECOND	[] [] [] []
		THIRD	[] [] [] []
Page : 2	Print	Prev Page	Next Page Save Return

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Item Name 48 FERR		Auto Rerun Condition (Absorbance)	
Auto Rerun SW		Absorbance Range	
On Off		Lower	On Off
Auto Rerun Range (Result)		Higher	On Off
On Off On Off		Prozone Range	
Lower Higher		On Off	
Serum	#		
Urine	#		
Plasma	#		
CSF	#		
Dialysis	#		
Other	#		
Bottle Size (ml)			
24 Items	36 Items		
Reagent1 60	Reagent 1 0		
Reagent2 R1 15.0	Reagent2 R1 0		
Reagent2 R2 8.0	Reagent2 R2 0		
Page : 3	Print	Prev Page	Save Return

Item Name 48 FERR						
Out-of-Range Table						
	NEAT		Pre Dilution		Auto-rerun Dilution (*10)	
	ABOVE	BELOW	ABOVE	BELOW	ABOVE	BELOW
	Cal#1	Cal#2-Cal#6	Cal#1	Cal#2-Cal#6	Cal#1	Cal#2-Cal#6
Serum	#	#	#	#	#	#
Urine	#	#	#	#	#	#
Plasma	#	#	#	#	#	#
CSF	#	#	#	#	#	#
Dialysis	#	#	#	#	#	#
Other	#	#	#	#	#	#
Page : 4			Prev Page	Save	Return	

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13 BIBLIOGRAPHIE

- Wick M, Pinggera W, Lehmann, P: Iron metabolism. Diagnosis and Therapy of Anemias, 3rd edition. Vienna/New York : Springer - Verlag, 1996.
- Albertini A, Arosio P, Chianaone E, Drysdale J (editors). Ferritins and isoferitins as biochemical markers. Amsterdam/New York/Oxford Elsevier, 1984.
- Finlayson NDC. Hereditary (primary) haemochromatosis. BMJ 1990; 301 : 350 - 351.
- Jacobs D, Kasten B, Demott W, and Wolfson W; Laboratory Test Handbook, 2nd Edition. Williams and Wilkins, 1990.
- Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication. WHO/DIL/LAB/99.1 Rev 2, 2002.
- Fairbanks, V.F. and Klee, C.G., Textbook of Clinical Chemistry, Tiez, N. W., (ed.), W.B. Saunder Co., Philadelphia, pp 1495-1588 (1986).
- Cook, J.D., Skikne, B.S., Lynch, S.R and Reusser, M.E., Blood, 68, 726 (1986).

Item Name 48 FERR						
Out-of-Range Table						
	NEAT		Pre Dilution		Auto-rerun Dilution(*10)	
	ABOVE	BELOW	ABOVE	BELOW	ABOVE	BELOW
	Cal#1	Cal#2-Cal#6	Cal#1	Cal#2-Cal#6	Cal#1	Cal#2-Cal#6
Serum	#	#	#	#	#	#
Urine	#	#	#	#	#	#
Plasma	#	#	#	#	#	#
CSF	#	#	#	#	#	#
Dialysis	#	#	#	#	#	#
Other	#	#	#	#	#	#
Page :	4		Prev Page	Save	Return	

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