

# Human IgA2 Kit for use on SPAPLUS®

For *in vitro* diagnostic use

Product Code: LK088.S

Product manufactured by:  
The Binding Site Group Ltd, 8 Calthorpe Road, Edgbaston, Birmingham, B15 1QT, U.K.  
www.bindingsite.co.uk  
Telephone : +44 (0)121 456 9500  
Fax : +44 (0)121 456 9749  
E-mail: info@bindingsite.co.uk

FDA (USA) Information:  
Analyte Name: Immunoglobulins A subclasses  
Complexity Category: Moderate

SPAPLUS® is a registered trademark of The Binding Site Group Limited (Birmingham, UK) in certain countries.



## 1 INTENDED USE

This kit is intended for quantifying human IgA subclass 2 (IgA2) in serum using the Binding Site SPAPLUS turbidimetric analyser. Measurement of this immunoglobulin aids in the diagnosis of abnormal protein metabolism and the body's lack of ability to resist infectious agents. The test result is to be used in conjunction with other clinical and laboratory findings.

## 2 SUMMARY AND EXPLANATION

Two classes of IgA have been identified in humans: IgA1, which accounts for 80-90% of total serum IgA, and IgA2 which is the major subclass in secretions such as milk. The two subclasses appear to be regulated independently. Antigenic sites on the IgA subclasses are responsible for the anaphylactic transfusion reactions experienced by some patients totally deficient in either IgA1 or IgA2. After repeated transfusion such patients may produce antibodies to these antigens. IgA subclasses are of further importance in that certain pathogenic microorganisms, including *Haemophilus influenzae*, are capable of enzymic cleavage of IgA1, leading to partial inactivation of this subclass. Recurrent sinopulmonary infections may therefore be related to deficiency of IgA2 which is resistant to these organisms (refs 1, 2 & 3).

## 3 PRINCIPLE

Evaluating the concentration of a soluble antigen by turbidimetry involves the addition of the test sample to a solution containing the appropriate antibody in a reaction vessel or cuvette. A beam of light is passed through the cuvette and, as the antigen-antibody reaction proceeds, the light passing through the cuvette is scattered increasingly as insoluble immune complexes are formed. Light scatter is monitored by measuring the decrease in intensity of the incident beam of light. The antibody in the cuvette is in excess so the amount of immune complex formed is proportional to the antigen concentration. A series of calibrators of known antigen concentration are assayed initially to produce a calibration curve of measured light scatter versus antigen concentration. Samples of unknown antigen concentration can then be assayed and the results read from the calibration curve. The sensitivity of turbidimetric assays can be increased by the use of particle enhancement (ref. 4). This entails linking the antibody to a suitably sized particle that increases the relative light-scattering signal of the antigen-antibody reaction.

## 4 REAGENTS

- 4.1 Human IgA2 Latex Reagent:** Consisting of monospecific sheep antibody coated onto polystyrene latex and is supplied in stabilised liquid form. Preservatives: 0.099% sodium azide, 0.1% E-amino-n-caproic acid (EACA), 0.01% Benzamide, 0.05% ProClin™ 300 and 1% Protease Inhibitor Cocktail 1 (Calbiochem Catalogue No. 539131).
- 4.2 Calibrator and Controls:** These consist of pooled human serum and are supplied in stabilised liquid form. They contain 0.099% sodium azide, 0.1% EACA and 0.01% benzamide as preservatives. The concentration of IgA2 given on the quality control certificate has been obtained by comparison with the DA470k international reference material.
- 4.3 Reaction Buffer:** Containing 0.099% sodium azide as a preservative.

ProClin™ is a trademark of Rohm and Haas Corp., Philadelphia, PA.

## 5 CAUTION

All donors of human serum supplied in this kit have been serum tested and found negative for hepatitis B surface antigen (HBsAg) and antibodies to human immunodeficiency virus (HIV1 and HIV2) and hepatitis C virus. The assays used were either approved by the FDA (USA) or cleared for *in vitro* diagnostic use in the EU (Directive 98/79/EC, Annex II); however, these tests cannot guarantee the absence of infective agents. **Proper handling and disposal methods should be established as for all potentially infective material, including (but not limited to) users wearing suitable gloves, protective equipment and clothing at all times.** Only personnel fully trained in such methods should be permitted to perform these procedures.

This product contains sodium azide and ProClin 300 and must be handled with caution. Do not ingest or allow contact with the skin (particularly broken skin or open wounds) or mucous membranes. If contact does occur wash with a large volume of water and seek medical advice. Explosive metal azides may be formed on prolonged contact of sodium azide with lead and copper plumbing; on disposal of reagent, flush with a large volume of water to prevent azide build up.

**This product should only be used by suitably trained personnel for the purposes stated in the Intended Use. Strict adherence to these instructions is essential at all times. Results are likely to be invalid if parameters other than those stated in these instructions are used.**

Reagents from different batch numbers of kits are **NOT** interchangeable. If large numbers of tests are performed care should be taken to ensure that all the reagents are from the same batch.

## 6 STORAGE AND STABILITY

The unopened kit should be stored at 2-8°C and can be used until the expiry date shown on the kit box label. DO NOT FREEZE. The Human IgA2 Latex Reagent, Reaction Buffer, Calibrators and Controls may be stored for up to three months after opening, provided that they are capped to avoid evaporation and kept at 2-8°C in a refrigerator. The Human IgA2 Latex Reagent and Reaction Buffer may be stored at 8-12°C uncapped, on the SPAPLUS analyser for up to 30 days, provided that the main power switch (located at the rear of the left hand panel) is left switched on.

## 7 SPECIMEN COLLECTION AND PREPARATION

Use fresh or deep frozen serum samples.

Blood samples should be collected by venepuncture, allowed to clot naturally and the serum separated as soon as possible to prevent haemolysis. The serum may be stored at 2-8°C for up to 48 hours prior to assay, or for prolonged storage kept at -20°C or below. Repeated freezing and thawing should be avoided. Microbially contaminated, haemolysed and lipaemic serum and samples containing particulate matter should not be used.

## 8 METHODOLOGY

### 8.1 Materials provided

- 8.1.1 1 x 50 Tests Human IgA2 Latex SPAPLUS
- 8.1.2 1 x Human IgA2 SPAPLUS Calibrator set 1-6 (6 x 1.0mL)
- 8.1.3 2 x 1.2mL Human IgA Subclass SPAPLUS High Control
- 8.1.4 2 x 1.2mL Human IgA Subclass SPAPLUS Low Control
- 8.1.5 1 x 50 Tests IgA2 Reaction Buffer SPAPLUS

### 8.2 Materials required but not provided

- 8.2.1 Equipment for collection and preparation of test samples e.g. sample tubes, centrifuge etc.
- 8.2.2 A fully operational and equipped SPAPLUS analyser.
- 8.2.3 Current analyser operating instructions: SPAPLUS Reference Guide, Insert Code FIN012.
- 8.2.4 Sample Diluent (99: Dil 1) Product Code: SN080.S

### 8.3 Reagent preparation

Before loading, gently mix by inversion ensuring no foam or bubbles are generated or remain on the surface as these may interfere with reagent aspiration.

### 8.4 Test procedure

The user should be familiar with the operation of the SPAPLUS analyser before attempting to carry out the test procedures. The analyser should be prepared for use according to the manufacturer's instructions and the assay protocol entered as described below.

For full details of analyser operation refer to the SPAPLUS Reference Guide, (FIN012) supplied with the analyser.

### 8.4.1 Test parameters

Assay parameters are entered into item number 12.

Item Name 12 IgA2		<b>CALIBRATION</b>		Auto Fill
<b>DATA INFORMATION</b>		Type	Logit 2	
Units	mg/L	Standard		
Decimals	1	1 #	4 #	
<b>ANALYSIS</b>		2 #	5 #	
Type	End	3 #	6 #	
Main W.Length 1	600	<b>NORMAL RANGE</b>		
Sub W.Length		LOW	HIGH	MALE
Method				FEMALE
<b>CORR.</b>		Serum	[ ]	[ ]
Y =	1 X + 0	Urine	[ ]	[ ]
		Plasma	[ ]	[ ]
		CSF	[ ]	[ ]
		Dialysis	[ ]	[ ]
		Other	[ ]	[ ]
Page : 1	Print	Hard Copy	Next Page	Save Return

Item Name 12 IgA2		<b>DATA PROCESS</b>		<b>ABSORBANCE LIMIT</b>
<b>ASPIRATION</b>		READ	START	END
KIND	<input type="radio"/> Single <input checked="" type="radio"/> Double	MAIN	53	54
VOLUME		SUB	35	36
<b>SAMPLE</b>	3	LOW	-3.0	
REAGENT1 VOL	145 µL	HIGH	3.0	
REAGENT2 VOL	80	<b>FACTOR</b>		
		Blank correction	1	<input type="radio"/> ON <input checked="" type="radio"/> OFF
		ENDPOINT LIMIT	2.0	CHECK POINT
		LINEAR CHECK (%)	0	LOW-3
				HIGH 3
Third mix	<input checked="" type="radio"/> OFF <input type="radio"/> ON	<b>DILUTION</b>		
R1 Blank	<input type="radio"/> Water <input checked="" type="radio"/> Blank	Diluent	<input checked="" type="radio"/> 99: Dil 1	<input type="radio"/> 100: Dil 2
		Pre Dilution Rate	10	
		Auto Rerun Dilution Rate High		
		Auto Rerun Dilution Rate Low		
<b>MONITOR</b>		<b>PROZONE CHECK</b>		
0 LEVEL SPAN	1	START	END	LIMIT (%)
SPAN	3.0	FIRST	[ ]	[ ]
		SECOND	[ ]	[ ]
		THIRD	[ ]	[ ]
			<input type="radio"/> Low	<input checked="" type="radio"/> High
			<input type="radio"/> Low	<input checked="" type="radio"/> High
Page : 2	Print	Hard Copy	Prev Page	Next Page Save Return

Item Name 12 IgA2		<b>Auto Rerun Condition (Absorbance)</b>	
<b>Auto Rerun SW</b>		<input type="radio"/> On <input checked="" type="radio"/> Off	
<b>Auto Rerun Range (Result)</b>		Absorbance Range	
<input checked="" type="radio"/> On <input type="radio"/> Off <input type="radio"/> On <input checked="" type="radio"/> Off		Lower <input checked="" type="radio"/> On <input type="radio"/> Off	
Lower Higher		Higher <input type="radio"/> On <input checked="" type="radio"/> Off	
Serum	Cal 1 #	Prozone Range <input type="radio"/> On <input checked="" type="radio"/> Off	
Urine	Cal 6 #		
Plasma			
CSF			
Dialysis			
Other			
<b>Bottle Size (ml)</b>			
24 Items	36 Items		
Reagent1 60	Reagent1		
Reagent2 R1 7.75	Reagent2 R1		
Reagent2 R2 4.5	Reagent2 R2		
Page : 3	Print	Prev Page	Save Return

NB. The calibrator (Standard #) values are found in the Quality Control Certificate (SIN179.QC). Calibrator values on **Page 1** should be entered in ascending order, i.e. the lowest value first. The analyser will automatically calculate and enter the correct measuring ranges on item pages 3 and 4 providing the **Autofill** button is pressed after typing the value for calibrator 6 on page 1. View Item parameter pages 3 and 4 to ensure correct value entry.

### 8.5 Measuring range

The approximate measuring range of the IgA2 assay when using the standard 1/10 sample dilution is 50-1250 mg/L.

## 9 QUALITY CONTROL

- 9.1 At least two levels of appropriate control material should be tested a minimum of once a day. In addition, controls should be tested after calibration, with each new lot of reagent and after specific maintenance or troubleshooting steps described in the SPAPLUS Operation Manual.
- 9.2 Quality control testing should be performed in accordance with regulatory requirements and each laboratory's standard procedure. Should a control measurement be out of range when assayed with a store curve the assay must be calibrated. If on recalibration the control values measured with the new curve are still out of range, the instrument and the assay parameters should be checked before repeating the assay. If problems persist, refer to the local technical support organisation.
- 9.3 The concentrations of the controls provided are stated on the accompanying QC certificate (SIN179.QC). Sample results obtained should only be accepted if the control results are within  $\pm 15\%$  of the concentration(s) stated.

## 10 LIMITATIONS

- 10.1 Turbidimetric assays are not suitable for measurement of highly lipaemic or haemolysed samples or samples containing high levels of circulating immune complexes (CICs) due to the unpredictable degree of non-specific scatter these sample types may generate. Unexpected results should be confirmed using an alternative assay method.
- 10.2 This assay has not been validated using paediatric samples.
- 10.3 Should a control measurement be out of range when assayed with a stored curve the assay must be recalibrated. If on recalibration the control values measured with the new curve are still out of range, the instrument and the assay parameters should be checked before repeating the assay. If problems persist, refer to supplier.
- 10.4 Diagnosis cannot be made and treatment must not be given on the basis of IgA2 measurements alone. Clinical history and other laboratory findings must be taken into account.
- 10.5 Variation in reagent temperature may affect results. Ensure that reagents are transferred directly from the refrigerator to the refrigerated reagent compartment of the analyser – do not allow to warm to room temperature.

## 11 EXPECTED VALUES

The ranges provided have been obtained from a limited number of samples and are intended for guidance purposes only. Wherever possible it is strongly recommended that local ranges are generated.

### Adult serum ranges

These ranges were obtained using this kit, by measuring the IgA2 concentration of sera taken from healthy adult UK blood donors. The reference interval was calculated using non-parametric statistics and represents the central 95% of the population.

	Number (n)	Mean (mg/L)	Median (mg/L)	95 Percentile Range (mg/L)
IgA2	120	392.5	321.9	68.9 – 1142.5

## 12 PERFORMANCE CHARACTERISTICS

### 12.1 Precision

A study was performed following CLSI *Evaluation of Precision Performance of Clinical Quantitative Measurement Methods; Approved Guideline* (CLSI Document EP5-A2). The study was performed over 21 working days, with two runs per day. One user assessed three different samples using three different reagent lots on three analysers.

	IgA2 Precision summary								
	Mean (mg/L)	Within run		Between run		Between day		Total	
		SD	CV %	SD	CV %	SD	CV %	SD	CV %
Serum 1	1121.47	31.14	2.8	16.43	1.5	68.36	6.1	76.89	6.9
Serum 2	684.56	15.46	2.3	10.92	1.6	25.95	3.8	32.11	4.7
Serum 3	79.93	3.37	4.2	3.01	3.8	6.00	7.5	7.51	9.4

### 12.2 Comparison

A correlation study was performed on 89 samples (26 normal serum and 64 clinical serum) using this kit on a SPAPLUS and an alternative commercially available IgA2 assay. The study demonstrated excellent agreement with the following Passing Bablok plot:

$$y = 0.99x - 1.16 \text{ (mg/L)} \quad (y = \text{SPAPLUS IgA2}; x = \text{alternative assay})$$

correlation coefficient  $r = 0.994$  (calculated by linear regression)

### 12.3 Limit of Blank and Limit of Detection

Based on CLSI document *EP17-A - Protocols for Determination of Limits of Detection and Limits Quantitation; Approved Guideline* the limit of detection represents the lowest measurable analyte level that can be distinguished from zero. This has been estimated at 1.87 mg/L (n = 60).

The limit of quantitation for this assay is defined as the lowest point of the calibration curve i.e. 4.9mg/L based upon measuring neat sample.

### 12.4 Linearity

A linearity study was performed based on NCCLS *Evaluation of the Linearity of Quantitative Measurement Procedures* (NCCLS Document EP6-A). One user assessed the linearity of a pool of high samples using one lot of reagent on one analyser. This gave a regression plot of  $y = 1.023x + 1.296$  (y = measured IgA2 concentration, x = theoretical concentration) over the range of 73.77- 1250mg/L using the analyser's 1/10 sample dilutions.

### 12.5 Interference

No significant assay interference by 1530 formazine turbidity units (FTU) of chyle, 200mg/L bilirubin, 4.83g/L haemoglobin has been demonstrated at the minimum sample dilution (1/1).

	Bilirubin	Hb	Chyle
Mean IgA2 (mg/L)	19.2	19.4	19.7
% interference	-0.93%	0.7%	1.0%

### 12.6 Antigen excess

No antigen excess was observed to a level of five times the top point of the assay; approximately 6.55g/L.

## 13 BIBLIOGRAPHY

- Protein Reference Unit Handbook of Clinical Immunochemistry (1990) Ed. A. Milford Ward, Publ. PRU Publications, Sheffield, 133-134.
- Mestecky, J. *et al* (1989). Comparative studies of the biological properties of human IgA subclasses. *Protides of the Biol. Fluids* 36, 179-182.
- Conley, M.E. *et al* (1983). Serum levels of IgA1 and IgA2 in children and in patients with IgA deficiency. *Mol. Immunol.* 20, 977-981.
- Tang LX, Showell P, Carr-Smith HD, Mead GP, Drew R and Bradwell AR (2000). Evaluation of F (ab')<sub>2</sub>-based latex-enhanced nephelometric reagents for free immunoglobulin light chains on the Behring Nephelometer™ II. *Clin. Chem* 46:6, Suppl. 2000, 705, pA181.