# Lipoprotein (a) Kit for use on SPAPLUS®

For *in-vitro* diagnostic use

Product code: LK098.S

The Binding Site Group Ltd, 8 Calthorpe Road, Edgbaston, Birmingham, B15 1QT, U.K.  
www.bindingsite.co.uk  
Telephone : +44 (0)121 456 9500  
Fax : +44 (0)121 456 9749  
E-mail: info@bindingsite.co.uk

In Europe and the USA, SPAPLUS® is a registered trademark of The Binding Site Group Ltd, Birmingham, UK.

FDA (USA) Information  
Analyte Name: Lipoprotein(a) (Lp(a))  
Complexity Cat.: Moderate

## 1 INTENDED USE

Immunoturbidimetric assay for the quantitative *in vitro* determination of Lipoprotein (a) (Lp (a)) in human serum or plasma. This product is suitable for use on the SPAPLUS analyser.

## 2 SUMMARY AND EXPLANATION

Lp (a) determination is intended for use in conjunction with clinical evaluation, patient risk assessment and other lipid tests to evaluate disorders of lipid metabolism and to assess coronary heart disease in specific populations.<sup>1</sup>

## 3 PRINCIPLE

Agglutination occurs due to an antigen-antibody reaction between Lp (a) in a sample and anti-Lp (a) antibody adsorbed to latex particles\*. This agglutination is detected as an absorbance change at 700 nm, proportional to the concentration of Lp (a) in the sample.  
\*Note: This product is licensed from Denka Seiken.

## 4 REAGENTS

- 4.1 **R1. Lipoprotein (a) Reaction Buffer** Containing glycine buffer 0.17M, Sodium chloride 1.08M, Sodium ethylenediamine tetra acetic acid disodium salt dihydrate 0.05M, Sodium azide ≤0.09% w/v.
- 4.2 **R2. Lipoprotein (a) Reagent** Suspension of latex particles coated with anti-Lp (a) antibodies, containing glycine 0.17M, Sodium chloride 0.1M, Sodium azide ≤0.09% w/v.
- 4.3 **Lipoprotein (a) Controls** Supplied at 2 levels, low and high. Lp (a) controls are supplied lyophilised.
- 4.4 **Lipoprotein (a) Calibrator 1-5** Material based on lyophilised human serum containing Lp (a). Calibration has been carried out and values have been assigned using an immunoturbidimetric method with reference material traceable to WHO Reference Material SRM2B.

## 5 CAUTION

Reagent: For *in vitro* diagnostic use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.

Health and Safety Data Sheets are available on request.

**The reagents must be used only for the purpose intended by suitably qualified laboratory personnel, under appropriate laboratory conditions.**

Calibrator Set and Controls: For *in vitro* diagnostic use only, do not pipette by mouth, exercise the normal precautions required for handling laboratory reagents.

This material has been tested for the HIV (Human Immunodeficiency Virus) Antibody, HBs Ag and HCV antibody and found to be non-reactive using FDA cleared methods. However, as no method can offer complete assurance as to the absence of infectious agents, this material should be handled as though capable of transmitting infectious disease.

Dispose of this material according to local regulations.

## 6 STORAGE AND STABILITY

### R1. Lipoprotein (a) Reaction Buffer

Buffer is ready for use and is stable up to the expiry date when stored at 2-8°C protected from light. The *Lipoprotein (a) Reaction Buffer* may be stored, uncapped, on the analyser for up to 30 days, provided that the main power switch (located at the rear of the left hand panel) is left switched on.

### R2. Lipoprotein (a) Reagent

The Lp (a) Latex Reagent is ready for use and stable up to the expiry date when stored at 2-8°C protected from light. Invert several times before use, avoiding the formation of foam. The *Lipoprotein (a) Reagent* may be stored, uncapped, on the analyser for up to 30 days, provided that the main power switch (located at the rear of the left hand panel) is left switched on.

### Lipoprotein (a) Calibrator Set 1-5 and Controls

Unreconstituted calibrator set and controls are stable up to the expiry date shown on the side of each individual bottle when stored at 2-8°C. Once reconstituted they are stable for 14 days at 2-8°C in the absence of bacterial contamination.

## 7 SPECIMEN COLLECTION AND PREPARATION

Collect serum using standard sampling tubes and plasma using tubes containing Li heparin or K EDTA<sup>-</sup>. Samples may be stored at 4°C for 14 days without significant decrease. For long term storage the samples should be stored at -20°C or -70°C.<sup>3</sup>

## 8 METHODOLOGY

### 8.1 Materials provided

- 8.1.1 1 x 100 tests *Lipoprotein (a) Reagent SPAPLUS*.
- 8.1.2 1 x 100 tests *Lipoprotein (a) Reaction Buffer SPAPLUS*.
- 8.1.3 1 x *Lipoprotein (a) SPAPLUS Calibrator 1-5* (5 x 1.0mL calibrator set for Lipoprotein (a)).
- 8.1.4 2 x 1.0mL *Lipoprotein (a) SPAPLUS Low Control*
- 8.1.5 2 x 1.0mL *Lipoprotein (a) SPAPLUS High Control*

### 8.2 Materials required but not provided

- 8.2.1 Equipment for collection and preparation of test samples e.g. sample tubes, centrifuge etc.
- 8.2.2 A fully operational and equipped SPAPLUS analyser.
- 8.2.3 Current analyser operating instructions: SPAPLUS Reference Guide, Insert Code FIN012.
- 8.2.4 Sample Diluent (99: Dil 1), Binding Site Product Code: SN080.S.

### 8.3 Calibrator Set, Controls and Reagent preparation

- 8.3.1 The calibrator set and controls are supplied in lyophilised form. Reconstitute each vial with exactly 1.0mL of distilled water and leave to stand for 30 minutes. Dissolve the contents completely by swirling or rotating. Ensure that no lyophilised material remains unreconstituted. Prior to use, mix the contents by inverting the vials. Note: The *Lipoprotein (a) SPAPLUS Calibrator* vials are not suitable for use with the SPAPLUS Calibrator Rack Direct Load (IK530).
- 8.3.2 Before loading the reagent, gently mix by inversion ensuring no foam or bubbles are generated or remain on the surface as these may interfere with reagent aspiration.

### 8.4 Test procedure

The user should be familiar with the operation of the SPAPLUS analyser before attempting to carry out the test procedures. The analyser should be prepared for use according to the manufacturer's instructions and the assay protocol entered as described below.

For full details of analyser operation refer to the SPAPLUS Reference Guide (FIN012) supplied with the analyser.

#### 8.4.1 Test parameters

Assay parameters are entered into item number 49.

Item Name 49 LPA		<b>DATA INFORMATION</b>		<b>CALIBRATION</b>	
Units	mg/dL	Type	Logit 2	Auto Fill	
Decimals	2	Standard		1 #	4 #
<b>ANALYSIS</b>		Type		2 #	5 #
Main W.Length 1		End		3 #	6 #
Sub W.Length		▼			
Method				<b>NORMAL RANGE</b>	
CORR.				MALE FEMALE	
Y =		SLOPE INTER		LOW HIGH	LOW HIGH
1 X + 0				Serum [ ] [ ] [ ] [ ] [ ] [ ]	
				Urine [ ] [ ] [ ] [ ] [ ] [ ]	
				Plasma [ ] [ ] [ ] [ ] [ ] [ ]	
				CSF [ ] [ ] [ ] [ ] [ ] [ ]	
				Dialysis [ ] [ ] [ ] [ ] [ ] [ ]	
				Other [ ] [ ] [ ] [ ] [ ] [ ]	
Page : 1		Print Hard Copy		Next Page Save Return	

Item Name 49 LPA		<b>DATA PROCESS</b>		<b>ASPIRATION</b>	
KIND		READ		START END	
Single Double		MAIN		53 54 LOW -3	
VOLUME		SUB		33 34 HIGH 3	
SAMPLE		REAGENT1 VOL		135 µL	
REAGENT2 VOL		65		<b>FACTOR</b>	
Third mix		Blank correction		Reaction Check	
Blank		Water - Blank		Blank correction ON OFF	
				ENDPOINT LIMIT 2 CHECK POINT	
				LINEAR CHECK (%) 0 LOW -3	
				HIGH 3	
		<b>DILUTION</b>		Diluent 99: Dil 1 100: Dil 2	
				Pre Dilution Rate	
				Auto Rerun Dilution Rate High 10	
				Auto Rerun Dilution Rate Low	
<b>MONITOR</b>		<b>PROZONE CHECK</b>			
0 LEVEL SPAN 1		FIRST		START END LIMIT (%) Min dOD [ 0 ]	
SPAN 3		SECOND		[ ] [ ] [ ] [ ] [ ] [ ]	
		THIRD		[ ] [ ] [ ] [ ] [ ] [ ]	
				Low High	
				Low High	
Page : 2		Print		Prev Page Next Page Save Return	

\*Automatically calculated

Item Name 49 LPA		<b>Auto Rerun SW</b>		<b>Auto Rerun Condition (Absorbance)</b>	
On Off		Lower Higher		Absorbance Range	
<b>Auto Rerun Range (Result)</b>		On Off		Lower Higher	
On Off On Off		On Off		On Off	
Serum # #		Prozone Range		On Off	
Urine # #					
Plasma # #					
CSF # #					
Dialysis # #					
Other # #					
<b>Bottle Size (ml)</b>					
24 Items		36 Items			
Reagent1 60		Reagent 1 0			
Reagent2 R1 14.5		Reagent2 R1 0			
Reagent2 R2 7.5		Reagent2 R2 0			
Page : 3		Print		Prev Page Next Page Save Return	

Item Name 49 LPA

**Out-of-Range Table**

	NEAT		Pre Dilution		Auto-rerun Dilution (*10)	
	ABOVE Cal#1	BELOW Cal#2-Cal#6	ABOVE Cal#1	BELOW Cal#2-Cal#6	ABOVE Cal#1	BELOW Cal#2-Cal#6
Serum	#	#	#	#	#	#
Urine	#	#	#	#	#	#
Plasma	#	#	#	#	#	#
CSF	#	#	#	#	#	#
Dialysis	#	#	#	#	#	#
Other	#	#	#	#	#	#

Page : 4 Prev Page Save Return

N.B. The calibrator (Standard #) values are found in the Quality Control certificate (SIN285.QC). Calibrator values on **Page 1** should be entered in ascending order, i.e. the lowest value first. **IMPORTANT:** the analyser will only update the calibrator values providing the **Auto Fill** button is pressed after typing the value for calibrator 5 on Page: 1. The assay utilises an extrapolated calibration curve therefore the Auto Rerun Range (Result) (Page 3) and Out-of-Range (Page 4) tables must be manually updated using values supplied in the Quality Control certificate (SIN285.QC).

**8.4.2 Control volumes**

Transfer 150µL of control fluid into a sample cup and place onto the rack.

**8.5 Measuring range**

The range of this assay is approximately 7 - 90mg/dL. In the event of a rerun, the upper limit of the assay range is extended to 180mg/dL. These values are dependent on the lot specific values of the calibrators in use.

SPAPLUS analyser dilution	Approximate measuring range (mg/dL)
1/1	7 - 90
1/10	70 - 180

**9 QUALITY CONTROL**

The controls provided should be included in all assay runs. The Lp (a) concentration is stated on the accompanying Quality Control certificate (SIN285.QC). Sample results obtained should only be accepted if the control results are within ±20% of the concentration(s) stated.


**10 LIMITATIONS**

- 10.1 These kits are not suitable for the measurement of lipaemic or haemolysed samples due to the unpredictable degree of non-specific scatter these sample types may generate. Unexpected results should be confirmed using an alternative assay method.
- 10.2 Customers are strongly advised to run all controls with every batch of samples being assayed. Should a control value be out of range against a stored curve, it is recommended that the control is re-assayed using the same calibration curve. If the control value is still out of range the curve should be recalibrated and the controls re-assayed. If control values are out of range against the new calibration curve check the instrument and parameters entered before repeating the assay. If problems persist, refer to the supplier.
- 10.3 Diagnosis cannot be made and treatment must not be initiated on the basis of Lp (a) measurements alone. Clinical history and other laboratory findings must also be taken into account.

**11 EXPECTED VALUES**

Adults	<30 mg/dL <sup>1,4,5</sup>
--------	----------------------------

The above reference range was established based on a sample of 96 Caucasian individuals comprising 49 males (age range 17-90 years; mean = 55 years) and 47 females (age range 13-84 years; mean = 55 years) resident in Northern Ireland. The population tested was an ambulatory population with no history of coronary disease. Results showed a mean Lp (a) value of 18.5mg/dL for males and 20.6mg/dL for females. Reference ranges have not been established for this assay for different ethnic populations or disease states.

Lp (a) concentrations have been shown to be genetically determined and to vary with ethnic populations. One study carried out in the United States showed that mean plasma levels of Lp (a) were approximately twice as high in African people or people of African descent compared to levels in Caucasians<sup>4</sup>. Also, the distribution of Lp (a) is less skewed in African people or people of African descent than in Caucasians<sup>4</sup>. Other studies have also shown no difference in Lp (a) levels between men (mean = 14mg/dL) and women (mean = 15mg/dL)<sup>5</sup>. Levels of Lp (a) have been shown not to differ significantly between pre- and post-menopausal Caucasian women<sup>5</sup>. It is therefore recommended that each laboratory establish its own reference range to reflect the age, race, sex, diet and geographical location of the population.

**12 PERFORMANCE CHARACTERISTICS**

**12.1 Precision**

	n	Mean (mg/dL)	Intra assay		Total precision	
			SD	CV %	SD	CV %
			Serum 1	80	20.03	0.57
Serum 2	80	46.82	0.63	1.3	1.87	4.0
Serum 3	80	70.64	0.61	0.9	2.36	3.3

**12.2 Comparison**

This method (Y) was compared with another commercially available method (X) and the following Passing & Bablok equation obtained:

$$y = 0.93x + 1.51$$

and a correlation coefficient of  $r = 0.996$

57 serum samples were analysed spanning the range 6.85 to 91.8 mg/dL.

**12.3 Analytical sensitivity**

The limit of Quantitation (LoQ), the limit of Detection (LoD) and the limit of Blank (LoB) were determined consistent with CLSI guidelines EP17-A2. LoQ is the smallest concentration that can be detected reliably. LoD is the smallest concentration that can be detected to determine the presence or absence of Lp (a). LoB is the highest concentration that is likely to be observed in a blank sample.

	Lp (a) concentration mg/dL
Limit of Blank	0.11
Limit of Detection	0.62
Limit of Quantitation	7.0

**12.4 Linearity**

The method is linear up to a Lp (a) concentration of 87.4 mg/dL.

This gave a regression plot of  $y = 0.97x + 0.71$ ;  $r = 0.998$  ( $y =$  measured Lp (a) concentration,  $x =$  theoretical concentration).

**12.5 Interference**

The analytes below were tested up to the following levels and were found not to interfere: Triglyceride 2000mg/dL, Haemoglobin 1000mg/dL, Free Bilirubin 60mg/dL, Conjugated Bilirubin 60mg/dL, Intralipid 1500mg/dL

**12.6 Antigen excess**

Antigen excess effects are not noted until levels approach 500 mg/dL.

**13 BIBLIOGRAPHY**

1. Lothar Thomas. Clinical Laboratory Diagnostics. 1st Edition. TH Books.
2. Tietz NW Clinical Guide to Laboratory Tests. Philadelphia, Pa: WB Sanders, 1995: 442 - 444.
3. Kronenberg, Lobentanz, Konig, Utermann and Dieplinger; (1994), Journal of Lipid Research 35:1318-1328.
4. Marcovina, S.M., Albers, J.J., Wijsman, E., Zhang, Z., Chapman, N.H. and Kennedy, H. (1996). Differences in Lp (a) concentrations and apo(a) polymorphs between black and white Americans. Journal of Lipid Research 37: 2569 - 2585
5. Jenner, J.L., Ordovas, J.m., Lamon-Fava, S. et al. (1993). Effects of Age, Sex and Menopausal Status on Plasma Lipoprotein (a) Levels. The Framingham Offspring Study. Circulation 87: 1135 - 1141.
6. Marcovina, S.M., Albers, J.J., Scanu, A.M. et al. (2000) Use of a Reference Material Proposed by the International Federation of Clinical Chemistry and Laboratory Medicine to Evaluate Analytical Methods for the Determination of Plasma Lipoprotein (a). Clinical Chemistry 46: (12) 1956 - 1967.

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## Zur *in-vitro* Diagnostik use

Bestell-Nr: LK098.S

The Binding Site Group Ltd., 8 Calthorpe Road, Edgbaston, Birmingham B15 1QT, UK.  
www.bindingsite.co.uk

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### 13 REFERENZEN

- Lothar Thomas. Clinical Laboratory Diagnostics. 1st Edition. TH Books.
- Tietz NW Clinical Guide to Laboratory Tests. Philadelphia, Pa: WB Sanders, 1995: 442 - 444.
- Kronenberg, Lobentanz, König, Utermann and Dieplinger; (1994), Journal of Lipid Research 35:1318-1328.
- Marcovina, S.M., Albers, J.J., Wijsman, E., Zhang, Z., Chapman, N.H. and Kennedy, H. (1996). Differences in Lp (a) concentrations and apo(a) polymorphs between black and white Americans. Journal of Lipid Research 37: 2569 – 2585
- Jenner, J.L., Ordovas, J.M., Lamon-Fava, S. et al. (1993). Effects of Age, Sex and Menopausal Status on Plasma Lipoprotein (a) Levels. The Framingham Offspring Study. Circulation 87: 1135 – 1141.
- Marcovina, S.M., Albers, J.J., Scanu, A.M. et al. (2000) Use of a Reference Material Proposed by the International Federation of Clinical Chemistry and Laboratory Medicine to Evaluate Analytical Methods for the Determination of Plasma Lipoprotein (a). Clinical Chemistry 46: (12) 1956 – 1967.

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The Binding Site Group Ltd, 8 Calthorpe Road, Edgbaston, Birmingham, B15 1QT, U.K.  
www.bindingsite.co.uk

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## 13 BIBLIOGRAPHIE

- Lothar Thomas. Clinical Laboratory Diagnostics. 1st Edition. TH Books.
- Tietz NW Clinical Guide to Laboratory Tests. Philadelphia, Pa: WB Sanders, 1995: 442 - 444.
- Kronenberg, Lobentanz, König, Utermann and Dieplinger; (1994), Journal of Lipid Research 35:1318-1328.
- Marcovina, S.M., Albers, J.J., Wijsman, E., Zhang, Z., Chapman, N.H. and Kennedy, H. (1996). Differences in Lp (a) concentrations and apo(a) polymorphs between black and white Americans. Journal of Lipid Research 37: 2569 - 2585
- Jenner, J.L., Ordovas, J.M., Lamon-Fava, S. et al. (1993). Effects of Age, Sex and Menopausal Status on Plasma Lipoprotein (a) Levels. The Framingham Offspring Study. Circulation 87: 1135 - 1141.
- Marcovina, S.M., Albers, J.J., Scanu, A.M. et al. (2000) Use of a Reference Material Proposed by the International Federation of Clinical Chemistry and Laboratory Medicine to Evaluate Analytical Methods for the Determination of Plasma Lipoprotein (a). Clinical Chemistry 46: (12) 1956 - 1967.

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- Lothar Thomas. Clinical Laboratory Diagnostics. 1st Edition. TH Books.
- Tietz NW Clinical Guide to Laboratory Tests. Philadelphia, Pa: WB Sanders, 1995: 442 - 444.
- Kronenberg, Lobentanz, Konig, Utermann and Dieplinger; (1994), Journal of Lipid Research 35:1318-1328.
- Marcovina, S.M., Albers, J.J., Wijsman, E., Zhang, Z., Chapman, N.H. and Kennedy, H. (1996). Differences in Lp (a) concentrations and apo(a) polymorphs between black and white Americans. Journal of Lipid Research 37: 2569 - 2585
- Jenner, J.L., Ordovas, J.m., Lamon-Fava, S. et al. (1993). Effects of Age, Sex and Menopausal Status on Plasma Lipoprotein (a) Levels. The Framingham Offspring Study. Circulation 87: 1135 - 1141.
- Marcovina, S.M., Albers, J.J., Scanu, A.M. et al. (2000) Use of a Reference Material Proposed by the International Federation of Clinical Chemistry and Laboratory Medicine to Evaluate Analytical Methods for the Determination of Plasma Lipoprotein (a). Clinical Chemistry 46: (12) 1956 - 1967.