# Rheumatoid Factor (RF) Kit for use on SPAPLUS®

For *in-vitro* diagnostic use

Product code: LK151.S

Manufactured in the UK  
The Binding Site Group Ltd, 8 Calthorpe Road, Edgbaston, Birmingham, B15 1QT, U.K.  
www.bindingsite.co.uk  
Telephone : +44 (0)121 456 9500  
Fax : +44 (0)121 456 9749  
E-mail: info@bindingsite.co.uk

In Europe and the USA, SPAPLUS® is a registered trademark of The Binding Site Group Ltd, Birmingham, UK.

## 1 INTENDED USE

A Rheumatoid Factor (RF) test system is a device intended for the quantitative *in vitro* determination of rheumatoid factor concentration in serum. This product is suitable for use on the SPAPLUS analyser.

## 2 SUMMARY AND EXPLANATION

This test system is intended to be used as an aid in the diagnosis of rheumatoid arthritis (RA). The agglutination reaction between antibody coated red blood cells and RA sera was first demonstrated by Waaler and Rose in 1940. The reaction has since been shown to be caused by certain factors in the RA sera. These rheumatoid factors are a heterogeneous group of high molecular weight auto-antibodies directed against the body's own immunoglobulins. They are produced by plasma cells present at sites of tissue injury. The initiating antigen is thought to be one or more viruses or viral antigens that persist in the joint tissues. Research has shown that both environmental and genetic factors can affect the production of rheumatoid factors with various biological properties. Rheumatoid factors have also been observed in the serum of patients with lupus erythematosus, hepatitis, liver cirrhosis, syphilis and various other conditions but the RF titre is much lower than in RA. Between 60 and 80% of patients with active RA possess this abnormal protein in their blood and joint fluid, and its detection is therefore of value in the diagnosis and monitoring of the disease.<sup>2</sup>

## 3 PRINCIPLE

Rheumatoid factors are antibodies directed against the Fc portion of IgG. The majority of rheumatoid factors are IgM antibodies, but may be IgG or IgA. Conditions giving rise to such factors include rheumatic conditions and chronic inflammatory processes. The RF latex reagent is a suspension of polystyrene latex particles of uniform size coated with human IgG. When serum containing rheumatoid factor is mixed with the RF latex reagent an increase in turbidity can be measured at 570nm. By constructing a standard curve from the absorbance of standards, the rheumatoid factor concentration can be determined<sup>1</sup>.

## 4 REAGENTS

- R1. RF Reaction Buffer** Containing Glycine buffer pH 8.3, Sodium Chloride, Sodium ethylenediamine tetra acetic acid disodium salt dihydrate, Bovine serum albumin, Sodium Azide 0.09% w/v.
- R2. RF Reagent** Containing Glycine Buffer pH 7.3, Sodium Chloride, Latex particle adsorbed human IgG, Sodium Azide 0.09% w/v.
- RF Controls** Supplied at 2 levels, Low and High. Target values and ranges are supplied in the Quality Control certificate. Supplied ready for use.
- RF Calibrator 1-5** Calibration has been carried out and values have been assigned using an immunoturbidimetric method standardised to the International Reference Preparation, WHO Standard 64/2. Supplied ready for use.

## 5 CAUTION

For *in vitro* diagnostic use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents. All reagents contain sodium azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention. Sodium azide reacts with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10% sodium hydroxide. Each human-based component in the kit has been found negative for the presence of HIV antibodies and hepatitis B surface antigen by FDA cleared methods. However, as no method can offer complete assurance, these reagents should be handled as if capable of transmitting infection. Health and Safety Data Sheets are available on request. Please dispose of all biological and chemical materials according to local guidelines. **The reagents must be used only for the purpose intended by suitably qualified laboratory personnel, under appropriate laboratory conditions.**

## 6 STORAGE AND STABILITY

- R1. RF Reaction Buffer**  
Supplied ready for use. Stable up to expiry date when stored at 2-8°C. The *RF Reaction Buffer* may be stored, uncapped, on the analyser for up to 35 days, provided that the main power switch (located at the rear of the left hand panel) is left switched on.
- R2. RF Reagent**  
Supplied ready for use. Stable up to expiry date when stored at 2-8°C. Invert several times before use, avoiding the formation of foam. The *RF Reagent* may be stored, uncapped, on the analyser for up to 35 days, provided that the main power switch (located at the rear of the left hand panel) is left switched on.
- RF Controls**  
OPENED: Store refrigerated (2-8°C). Protein control material is stable for 30 days at 2-8°C if kept capped in original container and free from contamination. Only the required amount of product should be removed. After use, any residual product should NOT BE RETURNED to the original vial.  
UNOPENED: Store refrigerated (2-8°C). Stable to expiration date printed on individual vials.

## RF Calibrator 1-5

OPENED: Calibrator material is stable for 28 days at 2-8°C if kept capped in original container and free from contamination. Only the required amount of product should be removed. After use, any residual product should NOT BE RETURNED to the original vial.  
UNOPENED: The RF Calibrators are supplied ready for use and are stable up to the expiry date when capped and stored at 2-8°C in the absence of contamination.

## 7 SPECIMEN COLLECTION AND PREPARATION

Fresh serum is recommended. The serum may be stored at 15-25°C for 24 hours or at 2-8°C for 3 days after collection. If the test cannot be carried out within this period of time, the serum should be frozen once at -15°C to -20°C for 4 weeks. Lipaemic or haemolytic samples and high levels of detergents may interfere in this assay.

## 8 METHODOLOGY

### 8.1 Materials provided

- 1 x 100 tests *RF Reagent SPAPLUS*.
- 1 x 100 tests *RF Reaction Buffer SPAPLUS*.
- 1 x *RF SPAPLUS Calibrator 1-5* (5 x 1.0mL calibrator set for RF).
- 1 x 1.0mL *RF SPAPLUS Low Control*.
- 1 x 1.0mL *RF SPAPLUS High Control*.

### 8.2 Materials required but not provided

- Equipment for collection and preparation of test samples e.g. sample tubes, centrifuge etc.
- A fully operational and equipped SPAPLUS analyser.
- Current analyser operating instructions: SPAPLUS Reference Guide, Insert Code FIN012.
- Sample Diluent (99: Dil 1) Product Code: SN080.S.

### 8.3 Reagent preparation

Before loading, gently mix by inversion ensuring no foam or bubbles are generated or remain on the surface as these may interfere with reagent aspiration.

### 8.4 Test procedure

The user should be familiar with the operation of the SPAPLUS analyser before attempting to carry out the test procedures. The analyser should be prepared for use according to the manufacturer's instructions and the assay protocol entered as described below.

For full details of analyser operation refer to the SPAPLUS Reference Guide (FIN012) supplied with the analyser.

### 8.4.1 Test parameters

Assay parameters are entered into item number 46

Item Name 46 RF		DATA INFORMATION		CALIBRATION	
Units	IU/mL	Type	Logit 2	Auto Fill	
Decimals	2	Standard			
ANALYSIS		1 #	4 #		
Type	End	2 #	5 #		
		3 #	6 #		
Main W.Length 1	570	NORM. RANGE			
Sub W.Length		MALE		FEMALE	
Method		LOW	HIGH	LOW	HIGH
CORR.		Serum	[ ]	[ ]	[ ]
Y =	1 X + 0	Urine	[ ]	[ ]	[ ]
		Plasma	[ ]	[ ]	[ ]
		CSF	[ ]	[ ]	[ ]
		Dialysis	[ ]	[ ]	[ ]
		Other	[ ]	[ ]	[ ]
Page : 1	Print	Hard Copy	Next Page	Save	Return

Item Name 46 RF		DATA PROCESS		ASPIRATION	
READ		START END		ABSORBANCE LIMIT	
KIND	Single	MAIN	45 46	LOW	-3
VOLUME	5	SUB	33 34	HIGH	2
SAMPLE REAGENT1 VOL	150	FACTOR			
REAGENT2 VOL	50	Blank correction	*	Reaction	OFF
		ENDPOINT LIMIT	2	CHECK POINT	
		LINEAR CHECK (%)	0	LOW	-3
				HIGH	3
Third mix	OFF	DILUTION			
Blank	Water - Blank	Diluent	99: Dil 1	100: Dil 2	
		Pre Dilution Rate			
		Auto Rerun Dilution Rate High	10		
		Auto Rerun Dilution Rate Low			
MONITOR		PROZONE CHECK			
LEVEL SPAN 1	3	START	END	LIMIT (%)	Min dOD [ 0 ]
		FIRST	[ ]	[ ]	
		SECOND	[ ]	[ ]	Low High
		THIRD	[ ]	[ ]	Low High
Page : 2	Print	Prev Page	Next Page	Save	Return

\*Automatically calculated

Item Name 46 RF		Auto Rerun SW		Auto Rerun Condition (Absorbance)	
Auto Rerun Range (Result)		Lower		Higher	
Serum #		On		Off	
Urine #		Off		On	
Plasma #		On		Off	
CSF #		Off		On	
Dialysis #		On		Off	
Other #		Off		On	
Bottle Size (ml)		Prozone Range		On	
24 Items	36 Items				
Reagent1 60	Reagent 1 0				
Reagent2 R1 16.0	Reagent2 R1 0				
Reagent2 R2 6.0	Reagent2 R2 0				
Page : 3	Print	Prev Page	Next Page	Save	Return

Item Name 46 RF						
<b>Out-of-Range Table</b>						
	NEAT		Pre Dilution		Auto-rerun Dilution(*10)	
	ABOVE	BELOW	ABOVE	BELOW	ABOVE	BELOW
	Cal#1	Cal#2-Cal#6	Cal#1	Cal#2-Cal#6	Cal#1	Cal#2-Cal#6
Serum	#	#	#	#	#	#
Urine	#	#	#	#	#	#
Plasma	#	#	#	#	#	#
CSF	#	#	#	#	#	#
Dialysis	#	#	#	#	#	#
Other	#	#	#	#	#	#
Page : 4			Prev Page	Save	Return	

N.B. The calibrator (Standard #) values are found in the Quality Control certificate (SIN282.QC). Calibrator values on Page 1 should be entered in ascending order, i.e. the lowest value first. IMPORTANT: the analyser will only update the calibrator values providing the **Auto Fill** button is pressed after typing the value for calibrator 5 on Page: 1. The assay utilises an extrapolated calibration curve therefore the Auto Rerun Range (Result) (Page 3) and Out-of-Range (Page 4) tables must be manually updated using values supplied in the Quality Control certificate (SIN282.QC).

#### 8.4.2 Control volumes

Transfer 150µL of control fluid into a sample cup and place onto the rack.

#### 8.5 Measuring range

The range of this assay is approximately 7 to 104 IU/mL. In the event of a rerun, the upper limit of the assay range is increased to 1040 IU/mL. These values are dependent on the lot specific value of the calibrator in use.

SPAPLUS analyser dilution	Approximate measuring range (IU/mL)
1/1	7 – 104
1/10	70 – 1040

### 9 QUALITY CONTROL

The controls provided should be included in all assay runs. The RF concentration is stated on the accompanying Quality Control certificate (SIN282.QC). Sample results obtained should only be accepted if the control results are within ±20% of the concentration(s) stated.

### 10 LIMITATIONS

- 10.1 These kits are not suitable for the measurement of samples containing circulating immune complexes (CICs) or for lipaemic or haemolysed samples due to the unpredictable degree of non-specific scatter these sample types may generate. Unexpected results should be confirmed using an alternative assay method.
- 10.2 Customers are strongly advised to run all controls with every batch of samples being assayed. Should a control value be out of range against a stored curve, it is recommended that the control is re-assayed using the same calibration curve. If the control value is still out of range the curve should be recalibrated and the controls re-assayed. If control values are out of range against the new calibration curve check the instrument and parameters entered before repeating the assay. If problems persist, refer to the supplier.
- 10.3 Diagnosis cannot be made and treatment must not be initiated on the basis of RF measurements alone. Clinical history and other laboratory findings must also be taken into account.

### 11 EXPECTED VALUES

The cut off value which constitutes the upper limit of normal is subject to dispute. However, values less than 10 IU/mL are considered to be within the normal range. The World Health Organisation has suggested a level approaching 12.5 IU/mL as being the upper limit of normal.<sup>1</sup> It is recommended that each laboratory should establish an expected range to reflect the geographical location and the age, sex and diet of the population.

### 12 PERFORMANCE CHARACTERISTICS

The following performance data was obtained using a SPAPLUS analyser.

#### 12.1 Precision

RF precision summary						
	n	Mean (IU/ml)	Intra assay		Total precision	
			SD	CV %	SD	CV %
Serum 1	80	15.28	0.58	3.8%	0.90	5.9%
Serum 2	80	43.77	1.36	3.1%	1.69	3.9%
Serum 3	80	76.7	1.87	2.4%	2.04	2.7%

#### 12.2 Comparison

This method (Y) was compared with another commercially available method (X) and the following Passing & Bablok equation obtained:

$$y = 0.98x - 0.39$$

and a correlation coefficient of  $r = 0.996$

68 patient samples were analysed spanning the range 7.2 to 99.5 IU/mL.

#### 12.3 Analytical sensitivity

The limit of Quantitation (LoQ), the limit of Detection (LoD) and the limit of Blank (LoB) were determined consistent with CLSI guidelines EP17-A2. LoQ is the smallest concentration that can be detected reliably. LoD is the smallest concentration that can be detected to determine the presence or absence of RF. LoB is the highest concentration that is likely to be observed in a blank sample.

#### SPAPLUS

	RF concentration IU/mL
Limit of Blank	2.3
Limit of Detection	3.9
Limit of Quantitation	7.0

#### 12.4 Linearity

The method is linear up to a RF concentration of 105.8 IU/mL.

This gave a regression plot of  $y = 1.03x - 0.58$   $r = 0.999$  ( $y$  = measured RF concentration,  $x$  = theoretical concentration).

#### 12.5 Interference

The analytes below were tested up to the following levels and were found not to interfere:  
 Haemoglobin 1000 mg/dL  
 Intralipid 2000 mg/dL  
 Free Bilirubin 60 mg/dL  
 Conjugated Bilirubin 60 mg/dL  
 Triglycerides 1500 mg/dL

#### 12.6 Antigen excess

Antigen excess effects are not noted until levels approach 390 IU/mL.

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Bestellnummer: LK151.S

In England hergestellt von:

The Binding Site Group Ltd., 8 Calthorpe Road, Edgbaston, Birmingham, B15 1QT, UK.  
www.bindingsite.co.uk

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