Human IgG CSF Kit for use on SPAPLUS

For *in-vitro* diagnostic use

Product Code: NK004.L.S

Product manufactured by:
The Binding Site Group Ltd., 8 Calthorpe Road, Edgbaston, Birmingham, B15 1QT, UK
www.bindingsite.co.uk
Telephone: +44 (0)121 456 9500
Fax: +44 (0)121 456 9749
E-mail: info@bindingsite.co.uk

FDA (USA) Information
Analyte Name: Immunoglobulins IgG
Complexity Cat.: Moderate

1 INTENDED USE

Human IgG CSF Kit for use on SPAPLUS is intended for the quantitative measurement of human IgG in cerebrospinal fluid (CSF) samples using the SPAPLUS analyser. Measurement of this immunoglobulin aids in the assessment of the body's lack of ability to resist infectious disease in conjunction with other clinical and laboratory findings.

2 SUMMARY AND EXPLANATION

Serum is the predominant source for proteins present in the CSF, the levels of which are regulated by the permeability of the blood-CSF barrier and CSF flow rate. An increase in CSF protein levels can be indicative of barrier dysfunction and/or local (intrathecal) synthesis of immunoglobulin (Ig) within the central nervous system (CNS)¹. These parameters can be evaluated by measurement of the serum and CSF concentrations of albumin, IgG, IgA and IgM.

As albumin in CSF originates exclusively from blood, the albumin CSF/serum ratio provides a measurement of barrier function. Calculation of CSF/serum ratios and comparison of the Ig ratios to the albumin CSF/serum value can differentiate between the serum-derived Ig and intrathecal Ig synthesis.

The assessment of barrier function, intrathecal synthesis and other variable CSF analytes can be useful in the diagnosis of a variety of CNS disorders^{1,2}.

3 PRINCIPLE

The determination of soluble antigen concentration by turbidimetric methods involves the reaction with specific antiserum to form insoluble complexes. When light is passed through the suspension formed a portion of the light is transmitted and focused onto a photodiode by an optical lens system. The amount of transmitted light is indirectly proportional to the specific protein concentration in the test sample. Concentrations are automatically calculated by reference to a calibration curve stored within the instrument.

4 REAGENTS

4.1 IgG CSF Antiserum: This is monospecific for IgG and is supplied in stabilised liquid form. It contains 0.099% sodium azide, 0.1% E-amino-n-caproic acid (EACA) and 0.01% benzamide as preservatives.

4.2 Calibrator and Controls: These consist of pooled human serum and are supplied in stabilised liquid form. They contain 0.099% sodium azide, 0.1% EACA and 0.01% benzamide as preservatives. The concentration of IgG given on the quality control certificate has been obtained by comparison with DA470k (formerly CRM470) International Reference Material.

4.3 Reaction Buffer: Containing 0.099% sodium azide as a preservative.

5 CAUTION

All donors of human serum supplied in this kit have been serum tested and found negative for hepatitis B surface antigen (HBsAg) and antibodies to human immunodeficiency virus (HIV1 and HIV2) and hepatitis C virus. The assays used were either cleared by the FDA (USA) or cleared for *in vitro* diagnostic use in the EU (Directive 98/79/EC, Annex II); however, these tests cannot guarantee the absence of infective agents. Proper handling and disposal methods should be established as for all potentially infective material, including (but not limited to) users wearing suitable protective equipment and clothing at all times. Only personnel fully trained in such methods should be permitted to perform these procedures.

WARNING: This product contains sodium azide and must be handled with caution; suitable gloves and other protective clothing should be worn at all times when handling this product. Do not ingest or allow contact with the skin (particularly broken skin or open wounds) or mucous membranes. If contact does occur wash with a large volume of water and seek urgent medical advice. Explosive metal azides may be formed on prolonged contact of sodium azide with lead and copper plumbing; on disposal of reagent, flush with a large volume of water to prevent azide build up.

This product should only be used by suitably trained personnel for the purposes stated in the Intended Use. Strict adherence to these instructions is essential at all times. Results are likely to be invalid if parameters other than those stated in these instructions are used.

Reagents from different batch numbers of kits are NOT interchangeable. If large numbers of tests are performed care should be taken to ensure that all the reagents are from the same batch.

6 STORAGE AND STABILITY

The unopened kit should be stored at 2-8°C and can be used until the expiry date shown on the kit box label. DO NOT FREEZE. The Human IgG CSF antiserum, Reaction Buffer, calibrators, and controls may be stored for up to two months after opening providing that they are capped to avoid evaporation and kept at 2-8°C in a refrigerator. The Human IgG

CSF antiserum and Reaction Buffer may be stored, uncapped, on the SPAPLUS analyser for up to 30 days, provided that the main power switch (located at the rear of the left hand panel) is left switched on.

7 SPECIMEN COLLECTION AND PREPARATION

Suitable samples are human CSF tested as fresh as possible. CSF samples may be stored at 2-8 °C for up to 7 days and can be kept at -20°C for up to 6 months.³ Samples must be centrifuged prior to testing.

8 METHODOLOGY

8.1 Materials provided

- 8.1.1 1 x 60 Tests Human IgG CSF Antiserum SPAPLUS
- 8.1.2 1 x Human IgG CSF SPAPLUS Calibrator set 1-6 (6 x 1.0mL)
- 8.1.3 2 x 1.5mL Combined CSF SPAPLUS High Control
- 8.1.4 2 x 1.5mL Combined CSF SPAPLUS Low Control
- 8.1.5 1 x 60 Tests IgG CSF Reaction Buffer SPAPLUS

8.2 Materials required but not provided

- 8.2.1 Equipment for collection and preparation of test samples e.g. sample tubes, centrifuge etc
- 8.2.2 A fully operational and equipped SPAPLUS analyser
- 8.2.3 Current analyser operating instructions: SPAPLUS Reference Guide, Insert Code FIN012
- 8.2.4 Sample Diluent (99: Dil 1) Product Code: SN080.S
- 8.2.5 2% Alkaline wash solution (working dilution)

8.3 Reagent preparation

Before loading, gently mix by inversion ensuring no foam or bubbles are generated or remain on the surface as these may interfere with reagent aspiration.

8.4 Test procedure

The user should be familiar with the operation of the SPAPLUS analyser before attempting to carry out the test procedures. The analyser should be prepared for use according to the manufacturer's instructions and the assay protocol entered as described below.

For full details of analyser operation refer to the SPAPLUS Reference Guide (FIN012) supplied with the analyser.

8.4.1 Test parameters

Assay parameters are entered into item number 30.

Item Name 30 IgG C		CALIBRATION	
DATA INFORMATION		Type Spline 1 ▼	
Units	mg/L	Standard	
Decimals	2	1 #	4 #
		2 #	5 #
		3 #	6 #
ANALYSIS			
Type	End ▼		
Main W.Length 1	340 ▼		
Sub W.Length	▼		
		NORMAL RANGE	
Method		LOW	MALE HIGH
			FEMALE LOW HIGH
CORR.		Serum	[] [] [] [] [] []
		Urine	[] [] [] [] [] []
		Plasma	[] [] [] [] [] []
		CSF	[] [] [] [] [] []
		Dialysis	[] [] [] [] [] []
		Other	[] [] [] [] [] []
Y =	SLOPE X + INTER		
	1 0		
Page : 1 Print Hard Copy		Next Page	Save Return

Item Name 30 IgG C		DATA PROCESS	
ASPIRATION		READ	ABSORBANCE LIMIT
KIND	Single Double	START	END
		MAIN	53 54
		SUB	30 31
		LOW	-3.0
		HIGH	#
SAMPLE		FACTOR	
VOLUME	30	Blank correction	* 0
REAGENT1 VOL	170 µL	ENDPOINT LIMIT	2.0
REAGENT2 VOL	20	CHECK POINT	
		LINEAR CHECK (%)	0
		LOW	-3
		HIGH	3
Third mix ● OFF ○ ON		DILUTION	
R1 Blank ● Water - Blank		Diluent	● 99: Dil 1 ○ 100: Dil 2
		Pre Dilution Rate	▼
		Auto Rerun Dilution Rate High	10
		Auto Rerun Dilution Rate Low	▼
MONITOR		PROZONE CHECK	
LEVEL SPAN 1		FIRST	START END LIMIT (%)
SPAN	3.0	SECOND	[] [] [] [] [] []
		THIRD	[] [] [] [] [] []
			○ Low ● High
			○ Low ● High
Page : 2 Print Hard Copy		Prev Page	Next Page Save Return

*Automatically calculated

Item Name 30 IgG C		Auto Rerun Condition (Absorbance)	
Auto Rerun SW		Absorbance Range	
● On ○ Off		Lower	○ On ● Off
Auto Rerun Range (Result)		Higher	○ On ● Off
○ On ● Off ● On ○ Off		Prozone Range	○ On ● Off
Lower Higher			
Serum	Cal 1 # Cal 6 #		
Urine			
Plasma			
CSF			
Dialysis			
Other			
Bottle Size (ml)			
24 Items	36 Items		
Reagent1	60	Reagent1	0
Reagent2 R1	10.9	Reagent2 R1	0
Reagent2 R2	1.9	Reagent2 R2	0
Page : 3 Print		Prev Page	Save Return

N.B. The calibrator (Standard #) and high absorbance limit (High #) values to be entered on Pages: 1 and 2 are found on the Quality Control Certificate (SIN213.QC). The calibrator values should be entered in ascending order, i.e. lowest value first.

12 PERFORMANCE CHARACTERISTICS

12.1 Precision

A study was performed following CLSI *Evaluation of Precision Performance of Clinical Quantitative Measurement Methods; Approved Guideline* (CLSI Document EP5-A2). The study was performed over 5 working days, with two runs per day. One user assessed three different samples using one reagent lot on one analyser.

	IgG CSF Precision Summary									
	Mean (mg/L)	Within run			Between run			Total		
		SD	CV %	SD	CV %	SD	CV %	SD	CV %	
Sample 1	125.533	1.3	1.0	6.2	4.9	0.0	0.0	6.3	5.0	
Sample 2	31.707	0.3	0.9	1.6	5.0	2.0	6.4	2.6	8.2	
Sample 3	7.447	0.1	1.6	0.4	5.5	0.2	3.3	0.5	6.6	

12.2 Comparison

A correlation study was performed on 96 CSF samples (35 normal CSF and 61 clinical CSF) using this kit on a SPAPLUS and an alternative commercially available IgG CSF assay. The study demonstrated the following Passing Bablok plot:

$$y = 1.05x - 0.27 \text{ (mg/L)} \quad (y = \text{SPAPLUS IgG CSF}; x = \text{alternative assay})$$

$$\text{correlation coefficient } r = 0.983 \quad (\text{calculated by linear regression})$$

12.3 Limit of Detection and Limit of Quantitation

Based on CLSI document EP17-A2 - *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline* the limit of detection represents the lowest measurable analyte level that can be distinguished from zero. This has been estimated at 0.69mg/L (n = 60).

The limit of quantitation has been calculated as 4.46mg/L (n=60) based on one lot at a 1/1 sample dilution.

12.4 Linearity

A linearity study was performed following CLSI (formerly NCCLS) *Evaluation of the Linearity of Quantitative Measurement Procedures* document EP6-A. One user assessed the linearity of a diluted high level pool of samples using one lot of reagent on one analyser. This gave a regression plot of $y = 0.9996x - 0.4893\text{mg/L}$ (y = measured IgG CSF concentration, x = theoretical concentration) over the range of 6.88 – 135.65mg/L.

12.5 Interference

No significant assay interference by 100mg/L bilirubin, 2.5g/L haemoglobin, 200mg/L acetaminophen and 600mg/L aspirin has been demonstrated at the minimum sample dilution (1/1).

	Bilirubin	Hb	Acetaminophen	Aspirin
Mean IgG (mg/L)	32.1	30.1	32.1	31.8
% interference	-2.59	-3.12	0.17	1.95

Bacterial interference has not been assessed (see section 10.7).

12.6 Antigen excess

No antigen excess was observed to a level of five times the top point of the assay; approximately 700mg/L.

13 BIBLIOGRAPHY

- Reiber H, Peter JB. Cerebrospinal fluid analysis: disease-related data patterns and evaluation programs. *J Neurol Sci* 2001;184:101-22.
- Regeniter A, Kuhle J, Mehling M, Möller H, Wurster U, Freidank H, Siede WH. A modern approach to CSF analysis: pathophysiology, clinical application, proof of concept and laboratory reporting. *Clin Neurol Neurosurg*. 2009 May;111(4):313-8.
- Wu AHB, ed. *Tietz Clinical guide to laboratory tests*, 4th ed. Philadelphia: WB Saunders Company; 2006: 606.
- Felgenhauer K. Laboratory diagnosis of neurological diseases. In Thomas L (Ed.) *Clinical laboratory diagnosis*, TH-books, Frankfurt/Main 1998; 1308-26.

IMPORTANT: the analyser will automatically calculate and enter the correct measuring ranges on Item pages 3 and 4 providing the **ENTER** button is pressed after typing the value for calibrator 6 on Page: 1. View Item parameter page 4 to ensure correct value entry.

8.4.2 Sample probe wash

To protect against any potential carry-over on the sample probe a sample probe wash must be programmed.

- 8.4.2.1 Click on **System** in the SPAPLUS software's Main Screen.
- 8.4.2.2 Select **Sys Para** from the drop down menu.
- 8.4.2.3 Select **Item** for the sample probe wash type and click **Save**. Click **OK** to save and **Exit** to return to the Main Screen.
- 8.4.2.4 Click on **System** in the SPAPLUS software's Main Screen.
- 8.4.2.5 Select **Sample Probe Wash** from the drop down menu.
- 8.4.2.6 Select **Item by item setting** and enter the parameters as below:

Sample Probe Wash (Item)			
Item 1	Item 2	Wash	
All clear	Line clear	W2	Update
			Exit

- 8.4.2.7 Click **Update** and **OK** to save the information and **Exit** to return to the Main Screen.

8.4.3 Calibration parameters

To protect against any potential carry-over on the sample probe from preceding assays six blank replicates must be programmed for the calibration curve.

Calibration Parameter							
CH ODR	ITEM#	Name	BLK ODR	Re CAL	BLK	STD-1	STD-2
1	30	IgG C			B1-6	S1-1	S2-1
<input type="checkbox"/>	Graph	030					
Order All						Update	Exit

8.4.4 Control volumes

For each control transfer 120µL into a sample cup and place on the rack.

8.4.5 Running controls and patient samples

Important: controls must be placed on the sample rack and not on the calibrator rack. If this is not followed the sample probe wash will not be activated.

- 8.4.5.1 Before placing the loaded sample rack onboard the analyser, fill a sample cup with 1.5mL of 2% alkaline wash solution (working dilution) and place into position W2 on the sample rack (W1 for albumin CSF). When sample testing is complete, discard the sample cup from the W2 position.
- 8.4.5.2 When ordering samples select the appropriate specimen type for each sample.

Note: After running CSF samples the 'sample type' in the **Order** screen will default to CSF. Ensure that 'specimen type' is changed as appropriate when manually ordering the next serum/urine sample.

8.5 Measuring range

The approximate measuring range of the assay is shown in the table below.

Specimen type	Analyser dilution	Approximate range mg/L
CSF	1/1	4.2 – 135
	1/10	42 - 1350

9 QUALITY CONTROL

The controls provided should be included in all assay runs. The IgG concentration is stated on the accompanying Quality Control Certificate (SIN213.QC). Sample results obtained should only be accepted if the control results are within ±15% of the concentration(s) stated.

10 LIMITATIONS

- Turbidimetric assays are not suitable for measurement of highly lipaemic or haemolysed samples or samples containing high levels of circulating immune complexes (CICs) due to the unpredictable degree of non-specific scatter these sample types may generate. Unexpected results should be confirmed using an alternative assay method.
- This assay has not been validated using paediatric samples.
- Should a control measurement be out of range when assayed with a stored curve the assay must be recalibrated. If on recalibration the control values measured with the new curve are still out of range, the instrument and the assay parameters should be checked before repeating the assay. If problems persist, refer to supplier.
- Diagnosis cannot be made and treatment must not be given on the basis of IgG measurements alone. Clinical history and other laboratory findings must be taken into account.
- Variation in reagent temperature may affect results. Ensure that reagents are transferred directly from the refrigerator to the refrigerated reagent compartment of the analyser – do not allow to warm to room temperature.
- Carry-over may occur in conditions where IgG levels are grossly elevated in sera e.g. with sera from multiple myeloma patients. Testing of such elevated samples should be isolated from CSF IgG testing.
- Bacterial interference has not been assessed. CSF samples should be as fresh as possible to limit bacterial growth and all samples must be centrifuged prior to testing (see section 7).

11 EXPECTED VALUES

The ranges provided have been obtained from a limited number of adult samples and are intended for guidance purposes only. Wherever possible it is strongly recommended that each facility should determine its own reference intervals since values may vary depending on the individual population studied.

Reference interval for IgG in CSF: <34mg/L (after conversion to DA470k).⁴

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Nur zur *in-vitro* Diagnostik

Bestell-Nr.: NK004.L.S

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The Binding Site Group Ltd., 8 Calthorpe Road, Edgbaston, Birmingham, B15 1QT, UK
www.bindingsite.co.uk

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Item Name 30 IgG C		CALIBRATION	
DATA INFORMATION			
Units	mg/L	Type	Spline 1 ▼
Decimals	2	Standard	
		1 #	4 #
		2 #	5 #
		3 #	6 #
ANALYSIS			
Type	End ▼		
Main W.Length 1	340 ▼		
Sub W.Length	▼		
NORMAL RANGE			
Method		LOW	HIGH
		MALE	FEMALE
		LOW	HIGH
		[]	[]
		[]	[]
		[]	[]
		[]	[]
		[]	[]
		[]	[]
		[]	[]
		[]	[]
CORR.			
Y =	1	X + 0	
SLOPE INTER			
Page : 1 Print Hard Copy Next Page Save Return			

Item Name 30 IgG C		DATA PROCESS	
READ			
ASPIRATION		START	END
KIND	Single	MAIN	53 54
	• Double	SUB	30 31
		LOW	-3.0
		HIGH	#
VOLUME			
SAMPLE	30	FACTOR	
REAGENT1 VOL	170 µL	Reaction Check	
REAGENT2 VOL	20	Blank correction * 0 ON • OFF	
		ENDPOINT LIMIT 2.0 CHECK POINT	
		LINEAR CHECK (%) 0 LOW 3	
		HIGH 3	
DILUTION			
Diluent		99: Dil 1	100: Dil 2
Pre Dilution Rate			
Auto Rerun Dilution Rate High	10		
Auto Rerun Dilution Rate Low			
PROZONE CHECK			
MONITOR			
0 LEVEL SPAN 1		START	END
SPAN	3.0	LIMIT (%)	
FIRST	[]	[]	[]
SECOND	[]	[]	[]
THIRD	[]	[]	[]
		Low	High
		Low	High
Page : 2 Print Hard Copy Prev Page Next Page Save Return			

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Item Name 30 IgG C		Auto Rerun SW		Auto Rerun Condition (Absorbance)	
Auto Rerun Range (Result)					
Absorbance Range					
Lower Higher • On • Off					
Prozone Range • On • Off					
Serum					
Urine					
Plasma					
CSF					
Dialysis					
Other					
Bottle Size (ml)					
24 Items 36 Items					
Reagent1	60	Reagent1	0		
Reagent2 R1	10.9	Reagent2 R1	0		
Reagent2 R2	1.9	Reagent2 R2	0		
Page : 3 Print Prev Page Save Return					

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13 REFERENZEN

1. Reiber H, Peter JB. Cerebrospinal fluid analysis: disease-related data patterns and evaluation programs. J Neurol Sci 2001;184:101-22.
2. Regeniter A, Kuhle J, Mehling M, Möller H, Wurster U, Freidank H, Siede WH. A modern approach to CSF analysis: pathophysiology, clinical application, proof of concept and laboratory reporting. Clin Neurol Neurosurg. 2009 May; 111(4): 313-8.
3. Wu AHB, ed. Tietz Clinical guide to laboratory tests, 4th ed. Philadelphia: WB Saunders Company; 2006: 606.
4. Felgenhauer K. Laboratory diagnosis of neurological diseases. In: Thomas L (Ed.) Clinical laboratory diagnosis, TH-books, Frankfurt/Main 1998; 1308-26.

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Item Name	30 IgG C
DATA INFORMATION	
Units	mg/L
Decimals	2
ANALYSIS	
Type	End ▼
Main W.Length 1	340 ▼
Sub W.Length	▼
Method	
CORR.	
SLOPE	INTER
Y =	1 X + 0
CALIBRATION	
Type	Spline 1 ▼
Standard	
1 #	4 #
2 #	5 #
3 #	6 #
NORMAL RANGE	
	MALE FEMALE
	LOW HIGH LOW HIGH
Serum	[] [] [] []
Urine	[] [] [] []
Plasma	[] [] [] []
CSF	[] [] [] []
Dialysis	[] [] [] []
Other	[] [] [] []
Page : 1	Print Hard Copy
	Next Page Save Return

Item Name	30 IgG C
DATA PROCESS	
READ	ASPIRATION
START END	MAIN 53 54 LOW -3.0
SUB 30 31	HIGH #3
VOLUME	SAMPLE
KIND <input type="radio"/> Single <input checked="" type="radio"/> Double	30
REAGENT1 VOL	170 µL
REAGENT2 VOL	20
FACTOR	Reaction Check
Blank correction	<input type="radio"/> ON <input checked="" type="radio"/> OFF
ENDPOINT LIMIT	* CHECK POINT
LINEAR CHECK (%) 0	LOW -3
	HIGH 3
DILUTION	
Diluent	<input checked="" type="radio"/> 99: Dil 1 <input type="radio"/> 100: Dil 2
Pre Dilution Rate	▼
Auto Rerun Dilution Rate High	10
Auto Rerun Dilution Rate Low	▼
MONITOR	PROZONE CHECK
0 LEVEL SPAN 1	START END LIMIT (%)
SPAN 3.0	FIRST [] [] [] <input type="radio"/> Low <input checked="" type="radio"/> High
	SECOND [] [] [] <input type="radio"/> Low <input checked="" type="radio"/> High
	THIRD [] [] [] <input type="radio"/> Low <input checked="" type="radio"/> High
Page : 2	Print Hard Copy
	Prev Page Next Page Save Return

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Item Name	30 IgG C
Auto Rerun SW	Auto Rerun Condition (Absorbance)
<input checked="" type="radio"/> On <input type="radio"/> Off	Absorbance Range
Auto Rerun Range (Result)	Lower <input type="radio"/> On <input checked="" type="radio"/> Off
<input type="radio"/> On <input checked="" type="radio"/> Off	Higher <input type="radio"/> On <input checked="" type="radio"/> Off
Serum Cal 1 #	Prozone Range <input type="radio"/> On <input checked="" type="radio"/> Off
Urine Cal 6 #	
Plasma	
CSF	
Dialysis	
Other	
Bottle Size (ml)	
24 Items	36 Items
Reagent1 60	Reagent1 0
Reagent2 R1 10.9	Reagent2 R1 0
Reagent2 R2 1.9	Reagent2 R2 0
Page : 3	Print
	Prev Page Save Return

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13 BIBLIOGRAPHIE

- Reiber H, Peter JB. Cerebrospinal fluid analysis: disease-related data patterns and evaluation programs. *J Neurol Sci* 2001;184:101-22.
- Regeniter A, Kuhle J, Mehling M, Möller H, Wurster U, Freidank H, Siede WH. A modern approach to CSF analysis: pathophysiology, clinical application, proof of concept and laboratory reporting. *Clin Neurol Neurosurg*. 2009 May; 111(4):313-8.
- Wu AHB, ed. Tietz Clinical guide to laboratory tests, 4th ed. Philadelphia: WB Saunders Company; 2006: 606.
- Felgenhauer K. Laboratory diagnosis of neurological diseases. In: Thomas L (Ed.) *Clinical laboratory diagnosis*, TH-books, Frankfurt/Main 1998; 1308-26.

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Item Name 30 IgG C		CALIBRATION	
DATA INFORMATION		Type Spline 1 ▼	
Units	mgL	Standard	
Decimals	2	1 #	4 #
ANALYSIS		2 #	5 #
Type	End ▼	3 #	6 #
Main W.Length 1	340 ▼	NORMAL RANGE	
Sub W.Length	▼		
Method		MALE	FEMALE
		LOW	HIGH
		LOW	HIGH
CORR.		Serum	[] [] [] [] [] []
		Urine	[] [] [] [] [] []
		Plasma	[] [] [] [] [] []
		CSF	[] [] [] [] [] []
		Dialysis	[] [] [] [] [] []
		Other	[] [] [] [] [] []
Y =	SLOPE INTER		
	1 X + 0		
Page : 1	Print	Hard Copy	Next Page Save Return

Item Name 30 IgG C		DATA PROCESS		ABSORBANCE LIMIT	
ASPIRATION		READ	START	END	LOW
KIND	Single Double	MAIN	53	54	3.0
		SUB	30	31	HIGH #
SAMPLE		VOLUME		FACTOR	
	30				Reaction Check
REAGENT1 VOL	170				Blank correction * ON OFF
REAGENT2 VOL	20				ENDPOINT LIMIT 2.0 CHECK POINT
					LINEAR CHECK (%) 0 LOW-3
					HIGH 3
Third mix	OFF ON	DILUTION			
R1 Blank	Water - Blank	Diluent		99: Dil 1	100: Dil 2
		Pre Dilution Rate			
		Auto Rerun Dilution Rate High	10		
		Auto Rerun Dilution Rate Low			
MONITOR		PROZONE CHECK			
0 LEVEL SPAN 1			START	END	LIMIT (%)
SPAN	3.0		FIRST		
			SECOND		
			THIRD		
Page : 2	Print	Hard Copy	Prev Page	Next Page	Save Return

Calculará automáticamente

Item Name 30 IgG C		Auto Rerun SW		Auto Rerun Condition (Absorbance)	
		On	Off	Absorbance Range	
Auto Rerun Range (Result)		Lower	Higher	Lower	On
		On	Off	Higher	On
Serum	Cal 1 #	Cal 6 #		Prozone Range	On
Urine					Off
Plasma					Off
CSF					Off
Dialysis					Off
Other					Off
Bottle Size (ml)					
24 Items		36 Items			
Reagent1	60	Reagent1	0		
Reagent2 R1	10.9	Reagent2 R1	0		
Reagent2 R2	1.9	Reagent2 R2	0		
Page : 3	Print		Prev Page	Save	Return

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