Apolipoprotein B (Apo B) Kit for use on SPAPLUS®

For *in-vitro* diagnostic use

Product code: NK086.S

The Binding Site Group Ltd, 8 Calthorpe Road, Edgbaston, Birmingham, B15 1QT, U.K.
www.bindingsite.co.uk
Telephone : +44 (0)121 456 9500
Fax : +44 (0)121 456 9749
E-mail: info@bindingsite.co.uk

In Europe and the USA, SPAPLUS® is a registered trademark of The Binding Site Group Ltd, Birmingham, UK.

FDA (USA) Information
Analyte Name: Apolipoprotein B
Complexity Cat.: Moderate

1 INTENDED USE

For the quantitative *in vitro* determination of apolipoprotein B (Apo B) in serum and plasma. This product is suitable for use on the SPAPLUS analyser.

2 SUMMARY AND EXPLANATION

Lipids are synthesised in the intestine or liver but must be transported to tissues and organs. However this is not possible without hydrophilic adaptation. Lipids are therefore transported by a series of micellar structures. These structures consist of an outer monolayer of protein (an apolipoprotein) and polar lipids (phospholipids and unesterified cholesterol) plus an inner core of neutral lipids (triglycerides and cholesterol esters). The apolipoproteins interact with a series of enzymes and tissue receptors and are therefore responsible for further metabolism and catabolism of the micelle. The B apolipoproteins are the main form of protein found in low density lipoproteins (LDL). Two forms of Apo B are found in humans. The most common form is B-100 (or large B) which constitutes the Apo B found in lipoproteins synthesised in the liver. The other form is B-48 (or small B), thought to be synthesised in the intestinal wall. Apo B is the main cholesterol carrying protein in the blood and is the ligand concerned with the uptake of cholesterol into cells by the LDL – receptor pathway. Apo B shows atherogenic signs and is thus useful for the evaluation of coronary risk. Studies have shown that there is an inverse relationship between Apo A-1 and coronary artery disease and a direct relationship with Apo B such that patients with CAD have generally reduced levels of Apo A-1 and increased levels of Apo B.^{1,2,3,4}

3 PRINCIPLE

This method is based on the reaction of a sample containing human Apo B and specific antiserum to form an insoluble complex which can be measured turbidimetrically at 340 nm. By constructing a standard curve from the absorbances of standards the concentration of Apo B can be determined.

4 REAGENTS

- 4.1 **R1. Apolipoprotein B Reaction Buffer** Containing polyethylene glycol <4%, Tris/HCl buffer 15mmol/L, Sodium chloride 106mmol/L.
- 4.2 **R2. Apolipoprotein B Antiserum** Anti-human Apo B.
- 4.3 **Apolipoprotein B Controls** Supplied at 2 levels, Medium and High. Target values and ranges are supplied in the Quality Control certificate. Supplied lyophilised.
- 4.4 **Apolipoprotein B Calibrator** Calibration has been carried out and values have been assigned using an immunoturbidimetric method with reference material standardised against WHO/IFCC Reference Standard SP3-07. Supplied lyophilised.

5 CAUTION

For *in vitro* diagnostic use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents. All reagents contain sodium azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention. Sodium azide reacts with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10% sodium hydroxide. Each human-based component in the kit has been found negative for the presence of HIV antibodies and hepatitis B surface antigen by FDA cleared methods. However, as no method can offer complete assurance, these reagents should be handled as if capable of transmitting infection. Health and Safety Data Sheets are available on request. Please dispose of all biological and chemical materials according to local guidelines. **The reagents must be used only for the purpose intended by suitably qualified laboratory personnel, under appropriate laboratory conditions.**

6 STORAGE AND STABILITY

- R1. Apolipoprotein B Reaction Buffer**
Supplied ready for use. Stable up to expiry date when stored at 2-8°C. The *Apolipoprotein B Reaction Buffer* may be stored, uncapped, on the analyser for up to 35 days, provided that the main power switch (located at the rear of the left hand panel) is left switched on.
- R2. Apolipoprotein B Antiserum**
Supplied ready for use. Stable up to expiry date when stored at 2-8°C. Invert several times before use, avoiding the formation of foam. The *Apolipoprotein B Reagent* may be stored, uncapped, on the analyser for up to 35 days, provided that the main power switch (located at the rear of the left hand panel) is left switched on.

Apolipoprotein B Controls

Unopened control are stable until the expiry date printed on the product label when stored between 2-8°C. Once reconstituted the components of the serum are stable for 7 days at 2-8°C, and 4 weeks at -20°C when frozen once.

Apolipoprotein B Calibrator

Unopened calibrator is stable until the expiry date printed on the product label when stored at 2-8°C. Once reconstituted Apolipoprotein B calibrators are stable for 7 days at 2-8°C, and 1 month at -20°C when frozen once.

7 SPECIMEN COLLECTION AND PREPARATION

APO B is stable for 1 day at 15-25°C or 8 days at 2-8°C and can be frozen once and stored for 2 months at -15°C to -25°C^{10,11}.

8 METHODOLOGY

8.1 Materials provided

- 8.1.1 1 x 100 tests *Apolipoprotein B Antiserum SPAPLUS*.
- 8.1.2 1 x 100 tests *Apolipoprotein B Reaction Buffer SPAPLUS*.
- 8.1.3 2 x *Apolipoprotein B SPAPLUS Calibrator* (2 x 1.0mL calibrator lyophilised).
- 8.1.4 2 x 1.0mL *Apolipoprotein B SPAPLUS Medium Control*
- 8.1.5 2 x 1.0mL *Apolipoprotein B SPAPLUS High Control*

8.2 Materials required but not provided

- 8.2.1 Equipment for collection and preparation of test samples e.g. sample tubes, centrifuge etc.
- 8.2.2 A fully operational and equipped SPAPLUS analyser.
- 8.2.3 Current analyser operating instructions: SPAPLUS Reference Guide, Insert Code FIN012.
- 8.2.4 Sample Diluent (99: Dil 1). Binding Site Product Code: SN080.S.

8.3 Calibrator, Control and Antiserum preparation

- 8.3.1 The calibrator and controls are supplied in lyophilised form. Open the vial carefully, avoiding any loss of the material and reconstitute with 1.0mL of distilled water. Replace the rubber stopper, close the vial and leave to stand for 30 minutes before use. Ensure that all traces of dry material are dissolved by swirling gently. Note: The *Apolipoprotein B SPAPLUS Calibrator* vials are not suitable for use with the SPAPLUS Calibrator Rack Direct Load (IK530).

A standard series must be manually produced following reconstitution of the Apolipoprotein B Calibrator fluid.

Standard	Dilution	Standard/diluent volumes (µl)
# 5	1/9.4	500 µL neat Calibrator fluid (reconstituted) + 4200 µL sample diluent
# 4	1/18.8	500 µL of Standard #5 + 500 µL diluent
# 3	1/37.6	500 µL of Standard #4 + 500 µL diluent
# 2	1/75.2	500 µL of Standard #3 + 500 µL diluent
# 1	1/150.4	500 µL of Standard #2 + 500 µL diluent

- 8.3.2 Antiserum: Before loading the reagent, gently mix by inversion ensuring no foam or bubbles are generated or remain on the surface as these may interfere with reagent aspiration.

8.4 Test procedure

The user should be familiar with the operation of the SPAPLUS analyser before attempting to carry out the test procedures. The analyser should be prepared for use according to the manufacturer's instructions and the assay protocol entered as described below.

For full details of analyser operation refer to the SPAPLUS Reference Guide (FIN012) supplied with the analyser.

8.4.1 Test parameters

Assay parameters are entered into item number 51

Item Name 51 APOB		CALIBRATION	
DATA INFORMATION		Type	Spline 1 ▼
Units	mg/dL	Standard	Auto Fill
Decimals	2	1 #	4 #
ANALYSIS		2 #	5 #
Type	End ▼	3 #	6 #
Main W.Length 1	340 ▼	NORMAL RANGE	
Sub W.Length	700 ▼	MALE HIGH FEMALE HIGH	
Method		Serum	[] [] [] []
CORR.		Urine	[] [] [] []
Y =	SLOPE INTER	Plasma	[] [] [] []
	1 X + 0	CSF	[] [] [] []
		Dialysis	[] [] [] []
		Other	[] [] [] []
Page : 1	Print Hard Copy	Next Page	Save Return

Item Name 51 APOB		DATA PROCESS	
ASPIRATION		READ	ABSORBANCE LIMIT
KIND	o Single • Double	MAIN 50 51	LOW -3
SAMPLE	VOLUME	SUB 30 31	HIGH 3
REAGENT1 VOL	200 µL	FACTOR	Reaction Check
REAGENT2 VOL	40	Blank correction *	o ON • OFF
		ENDPOINT LIMIT 2	CHECK POINT
		LINEAR CHECK (%) 0	LOW -3
Third mix	• OFF o ON		HIGH 3
Blank	• Water - Blank	DILUTION	
		Diluent	• 99: Dil 1 o 100: Dil 2
		Pre Dilution Rate	20 ▼
		Auto Rerun Dilution Rate High	100 ▼
		Auto Rerun Dilution Rate Low	▼
MONITOR		PROZONE CHECK	
0 LEVEL SPAN 1		START END LIMIT (%) Min dOD [0]	
SPAN 3		FIRST [] [] [] []	
		SECOND [] [] [] []	o Low • High
		THIRD [] [] [] []	o Low • High
Page : 2	Print	Prev Page	Next Page Save Return

*Automatically calculated

Item Name 51 APOB	
Auto Rerun SW	
<input type="radio"/> On <input type="radio"/> Off <input type="radio"/> On <input type="radio"/> Off	
Auto Rerun Range (Result)	
<input type="radio"/> On <input type="radio"/> Off <input type="radio"/> On <input type="radio"/> Off <input type="radio"/> Lower <input type="radio"/> Higher	
Serum #	#
Urine #	#
Plasma #	#
CSF #	#
Dialysis #	#
Other #	#
Bottle Size (ml)	
24 Items	36 Items
Reagent1 60	Reagent 1 0
Reagent2 R1 21.0	Reagent2 R1 0
Reagent2 R2 5.0	Reagent2 R2 0
Page : 3	Print Prev Page Save Return

Item Name 51 APOB						
Out-of-Range Table						
	NEAT	Pre Dilution (*20)		Auto-rerun Dilution (*100)		
	ABOVE	BELOW	ABOVE	BELOW	ABOVE	BELOW
	Cal#1	Cal#2-Cal#6	Cal#1	Cal#2-Cal#6	Cal#1	Cal#2-Cal#6
Serum #	#	#	#	#	#	#
Urine #	#	#	#	#	#	#
Plasma #	#	#	#	#	#	#
CSF #	#	#	#	#	#	#
Dialysis #	#	#	#	#	#	#
Other #	#	#	#	#	#	#
Page : 4	Prev Page		Save	Return		

N.B. The calibrator (Standard #) values are found in the Quality Control certificate (SIN287.QC). Calibrator values on Page 1 should be entered in ascending order, i.e. the lowest value first. **IMPORTANT:** the analyser will only update the calibrator values providing the **Auto Fill** button is pressed after typing the value for calibrator 5 on Page: 1. The assay utilises an extrapolated calibration curve therefore the Auto Rerun Range (Result) (Page 3) and Out-of-Range (Page 4) tables must be manually updated using values supplied in the Quality Control certificate (SIN287.QC).

8.4.2 Control volumes

Transfer 150µL of control fluid into a sample cup and place onto the rack.

8.5 Measuring range

The range of this assay is approximately 11.2 - 184 mg/dL. In the event of a rerun, the upper limit of the assay range is increased to approximately 920 mg/dL. These values are dependent on the lot specific value of the calibrator in use.

SPAPLUS analyser dilution	Approximate measuring range (mg/dL)
1/20	11.2 - 184
1/100	56 - 920

9 QUALITY CONTROL

The controls provided should be included in all assay runs. The Apolipoprotein B concentration is stated on the accompanying Quality Control certificate (SIN287.QC). Sample results obtained should only be accepted if the control results are within ±20% of the concentration(s) stated.

10 LIMITATIONS

- 10.1 These kits are not suitable for the measurement of samples containing circulating immune complexes (CICs) or for lipaemic or haemolysed samples due to the unpredictable degree of non-specific scatter these sample types may generate. Unexpected results should be confirmed using an alternative assay method.
- 10.2 Customers are strongly advised to run all controls with every batch of samples being assayed. Should a control value be out of range against a stored curve, it is recommended that the control is re-assayed using the same calibration curve. If the control value is still out of range the curve should be recalibrated and the controls re-assayed. If control values are out of range against the new calibration curve check the instrument and parameters entered before repeating the assay. If problems persist, refer to the supplier.
- 10.3 Diagnosis cannot be made and treatment must not be initiated on the basis of apolipoprotein B measurements alone. Clinical history and other laboratory findings must also be taken into account.

11 EXPECTED VALUES

NORMAL RANGE⁵

63 - 114 mg/dL.

It is recommended that each laboratory establish its own reference range to reflect the age, sex, diet and geographical location of the population.

12 PERFORMANCE CHARACTERISTICS

12.1 Precision

Apolipoprotein B precision summary						
	n	Mean (mg/dL)	Intra assay		Total precision	
			SD	CV %	SD	CV %
Serum 1	80	30.86	1.04	3.4%	2.22	7.2%
Serum 2	80	70.08	2.41	3.4%	3.26	4.7%
Serum 3	80	148.43	3.06	2.1%	4.17	2.8%

12.2 Comparison

This method (Y) was compared with another commercially available method (X) and the following linear regression equation obtained:

$$y = 1.05x - 7.15$$

and a correlation coefficient of $r = 0.990$

98 patient samples were analysed spanning the range 15.1 to 209.2 mg/dL.

12.3 Analytical sensitivity

The limit of Quantitation (LoQ), the limit of Detection (LoD) and the limit of Blank (LoB) were determined consistent with CLSI guidelines EP17-A2. LoQ is the smallest concentration that can be detected reliably. LoD is the smallest concentration that can be detected to determine the presence or absence of Apolipoprotein B. LoB is the highest concentration that is likely to be observed in a blank sample.

SPAPLUS

Apolipoprotein B concentration mg/dL	
Limit of Blank	2.08
Limit of Detection	3.81
Limit of Quantitation	9.00

12.4 Linearity

The method is linear up to an Apolipoprotein B concentration of 206.9 mg/dL.

This gave a regression plot of $y = 0.96x + 4.23$ $r = 0.998$ ($y =$ measured Apolipoprotein B concentration, $x =$ theoretical concentration).

12.5 Interference

The analytes below were tested up to the following levels and were found not to interfere:

- Triglyceride 2000 mg/dL
- Haemoglobin 1000 mg/dL
- Free Bilirubin 60 mg/dL
- Conjugated Bilirubin 60 mg/dL
- Intralipid 250 mg/dL

12.6 Antigen excess

Antigen excess effects are not noted until levels approach 3420 mg/dL.

13 BIBLIOGRAPHY

1. Bachorik, P.S., Kwiterovich, P.O. Clinica Chimica Acta 1988; 178: 1-34
2. Wang, X. L., Dudman, N.P.B. Clin. Chem. 1989; 35(10): 2082-2086
3. Wang, X. L., Dudman, N.P.B. Wilken, D.E.L., Clin Chem. 1989; 35 (6): 1000 - 1004.
4. Kukita H., Hiwada K., Kobubu T. Serum Apolipoprotein A-1, A-2 and B levels and their discriminative values in relatives of patients with coronary artery disease. Atheroscler 51: 261, 1984.
5. Provisional normal values recalculated on the basis of the CDC values C. Fruchart, J.-C. (1986), Ann. Biol. Clin. 44:116.
6. Adolphson, J.L., Albers J.J., Journal of Lipid Research 1989; 30: 599-806
7. Labour, C., Shepherd, J., Rosseneu, M. Clin. Chem. 1990; 36(4): 591-597
8. Marcynova, S.M, et al (1992), WHO/IFCC Meeting On Standardisation of Apolipoproteins, May 1992, Nice, France.
9. Albers, J.J, et al (1992) Clin Chem. 38:58.
10. Evans K, Mitcheson J, Laker M. Effect of Storage at 4°C and -20°C on Lipid, Lipoprotein, and Apolipoprotein Concentrations. Clin Chem 1995; 41: 392-396.
11. Use of anticoagulants in Diagnostic laboratory Investigations. WHO Publication. WHO/DIL/LAB 99.1/Rev. 2 Jan. 2002.

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13 REFERENZEN

- Bachorik, P.S., Kwiterovich, P.O. Clinica Chimica Acta 1988; 178: 1-34
- Wang, X.L., Dudman, N.P.B. Clin. Chem. 1989; 35(10): 2082-2086
- Wang, X. L., Dudman, N.P.B. Wiliken, D.E.L., Clin Chem. 1989; 35 (6): 1000 - 1004.
- Kukita H., Hiwada K., Kobubu T. Serum Apolipoprotein A-1, A-2 and B levels and their discriminative values in relatives of patients with coronary artery disease. Atheroscler 51: 261, 1984.
- Provisional normal values recalculated on the basis of the CDC values C. Fruchart, J-C. (1986), Ann. Biol. Clin. 44:116.
- Adolphson, J.L., Albers J.J., Journal of Lipid Research 1989; 30: 599-806
- Labour, C., Shepherd. J., Rosseneu, M. Clin. Chem. 1990; 36(4): 591-597
- Marcinova, S.M, et al (1992), WHO/IFCC Meeting On Standardisation of Apolipoproteins, May 1992, Nice, France.
- Albers, J.J, et al (1992) Clin Chem. 38:58.
- Evans K, Mitcheson J, Laker M. Effect of Storage at 4°C and -20°C on Lipid, Lipoprotein, and Apolipoprotein Concentrations. Clin Chem 1995; 41: 392-396.
- Use of anticoagulants in Diagnostic laboratory Investigations. WHO Publication. WHO/DIL/LAB 99.1/Rev. 2 Jan. 2002.

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The Binding Site Group Ltd, 8 Calthorpe Road, Edgbaston, Birmingham, B15 1QT, U.K.
www.bindingsite.co.uk

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Item Name 51 APOB		
DATA INFORMATION	CALIBRATION	
Units mg/dL	Type Spline 1 ▼	Auto Fill
Decimals 2	Standard	
	1 # 4 #	
ANALYSIS	2 # 5 #	
Type End ▼	3 # 6 #	
Main W.Length 1 340 ▼		
Sub W.Length 700 ▼		
	NORMAL RANGE	
Method	LOW HIGH MALE LOW HIGH FEMALE	
CORR.	Serum [] [] [] [] [] []	
	Urine [] [] [] [] [] []	
	Plasma [] [] [] [] [] []	
	CSF [] [] [] [] [] []	
	Dialysis [] [] [] [] [] []	
	Other [] [] [] [] [] []	
Y =	1 X + 0	
Page : 1	Print Hard Copy	Next Page Save Return

Item Name 51 APOB		
DATA PROCESS	READ	ABSORBANCE LIMIT
ASPIRATION	START END	
KIND ○ Single ● Double	MAIN 50 51	LOW -3
	SUB 30 31	HIGH 3
SAMPLE	VOLUME	
REAGENT1 VOL 200 µL		
REAGENT2 VOL 40 µL		
FACTOR	Reaction Check	
Blank correction *	○ ON ● OFF	
ENDPOINT LIMIT 2	CHECK POINT	
LINEAR CHECK (%) 0	LOW HIGH	-3 3
Third mix ● OFF ○ ON		
Blank ● Water - Blank		
	DILUTION	
Diluent	● 99: Dil 1 ○ 100: Dil 2	
Pre Dilution Rate 20		
Auto Rerun Dilution Rate High 100		
Auto Rerun Dilution Rate Low		
MONITOR	PROZONE CHECK	
0 LEVEL SPAN 1	START END LIMIT (%) Min dOD [0]	
SPAN 3	FIRST [] [] [] [] [] []	
	SECOND [] [] [] [] [] []	○ Low ● High
	THIRD [] [] [] [] [] []	○ Low ● High
Page : 2	Print	Prev Page Next Page Save Return

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Item Name 51 APOB

Auto Rerun SW
 On Off

Auto Rerun Range (Result)
 On Off

Lower Higher
 Lower Higher

Serum	#	#
Urine	#	#
Plasma	#	#
CSF	#	#
Dialysis	#	#
Other	#	#

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13 BIBLIOGRAPHIE

- Bachorik, P.S., Kwiterovich, P.O. Clinica Chimica Acta 1988; 178: 1-34
- Wang, X.L., Dudman, N.P.B. Clin. Chem. 1989; 35(10): 2082-2086
- Wang, X. L., Dudman, N.P.B. Wiliken, D.E.L., Clin Chem. 1989; 35 (6): 1000 - 1004.
- Kukita H., Hiwada K., Kobubu T. Serum Apolipoprotein A-1, A-2 and B levels and their discriminative values in relatives of patients with coronary artery disease. Atheroscler 51: 261, 1984.
- Provisional normal values recalculated on the basis of the CDC values C. Fruchart, J.-C. (1986), Ann. Biol. Clin. 44:116.
- Adolphson, J.L., Albers J.J., Journal of Lipid Research 1989; 30: 599-806
- Labour, C., Shepherd, J., Rosseneu, M. Clin. Chem. 1990; 36(4): 591-597
- Marcinova, S.M, et al (1992), WHO/FCC Meeting On Standardisation of Apolipoproteins, May 1992, Nice, France.
- Albers, J.J, et al (1992) Clin Chem. 38:58.
- Evans K, Mitcheson J, Laker M. Effect of Storage at 4°C and -20°C on Lipid, Lipoprotein, and Apolipoprotein Concentrations. Clin Chem 1995; 41: 392-396.
- Use of anticoagulants in Diagnostic laboratory Investigations. WHO Publication. WHO/DIL/LAB 99.1/Rev. 2

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Item Name 51 APOB		
DATA PROCESS		
ASPIRATION		
KIND	Single	Double
VOLUME		
SAMPLE 12		
REAGENT1 VOL 200 µL		
REAGENT2 VOL 40		
FACTOR		
Blank correction *		
ENDPOINT LIMIT 2		
LINEAR CHECK (%) 0		
Reaction Check		
Blank	OFF	ON
Water - Blank		
DILUTION		
Diluent	99: Dil 1	100: Dil 2
Pre Dilution Rate	20	
Auto Rerun Dilution Rate High	100	
Auto Rerun Dilution Rate Low		
MONITOR		
0 LEVEL SPAN 1		
SPAN 3		
PROZONE CHECK		
START	END	LIMIT (%)
FIRST [] [] []		Min dOD [0]
SECOND [] [] []		Low High
THIRD [] [] []		Low High
Page : 2	Print	Prev Page
	Next Page	Save
		Return

* Calculado automáticamente

Item Name 51 APOB

Auto Rerun SW
 On Off

Auto Rerun Range (Result)
 On Off On Off

Lower Higher
 Serum # #
 Urine # #
 Plasma # #
 CSF # #
 Dialysis # #
 Other # #

Auto Rerun Condition (Absorbance)
Absorbance Range
 Lower On Off
 Higher On Off

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13 BIBLIOGRAFIA

- Bachorik, P.S., Kwiterovich, P.O. Clinica Química Acta 1988; 178: 1-34
- Wang, X.L., Dudman, N.P.B. Clin. Chem. 1989; 35(10): 2082-2086
- Wang, X. L., Dudman, N.P.B. Wilken, D.E.L., Clin Chem. 1989; 35 (6): 1000 - 1004.
- Kukita H., Hiwada K., Kobubu T. Serum Apolipoprotein A-1, A-2 and B levels and their discriminative values in relatives of patients with coronary artery disease. Atheroscler 51: 261, 1984.
- Provisional normal values recalculated on the basis of the CDC values C. Fruchart, J-C. (1986), Ann. Biol. Clin. 44:116.
- Adolphson, J.L., Albers J.J., Journal of Lipid Research 1989; 30: 599-806
- Labour, C., Shepherd, J., Rosseneu, M. Clin. Chem. 1990; 36(4): 591-597
- Marcinova, S.M, et al (1992), WHO/IFCC Meeting On Standardisation of Apolipoproteins, May 1992, Nice, France.
- Albers, J.J, et al (1992) Clin Chem. 38:58.
- Evans K, Mitcheson J, Laker M. Effect of Storage at 4°C and -20°C on Lipid, Lipoprotein, and Apolipoprotein Concentrations. Clin Chem 1995; 41: 392-396.
- Use of anticoagulants in Diagnostic laboratory Investigations. WHO Publication. WHO/DIL/LAB 99.1/Rev. 2