

DHEA-S ELISA

Enzyme immunoassay for the quantitative determination of
DHEA-S in human serum or plasma.

REF **RE52181**

 **1 x 96**

   **2-8°C**

EU: **IVD**  U.S.: *For in-vitro
diagnostic use only.*



I B L I N T E R N A T I O N A L G M B H

Flughafenstrasse 52a
D-22335 Hamburg, Germany

Phone: +49 (0)40-53 28 91-0
Fax: +49 (0)40-53 28 91-11

IBL@IBL-International.com
www.IBL-International.com

1. INTENDED USE

Enzyme immunoassay for the *in-vitro diagnostic* quantitative determination of dehydroepiandrosterone sulfate (DHEA-S) in human serum or plasma.

2. SUMMARY AND EXPLANATION

Dehydroepiandrosterone sulfate (DHEA-S), a C-19 steroid hormone, is the most abundant adrenal androgen in the circulation. Most of the circulating DHEA-S originates from either direct adrenal secretion or by peripheral sulfation of DHEA secreted by the adrenal cortex. Because it is sulfated, it has a long half-life and hence lacks a circadian rhythm. DHEA-S does not circulate bound to specific proteins. It circulates at much higher levels than other androgens or related steroids. Levels of DHEA-S increase from about the seventh year of life, peak in the third decade and decrease gradually thereafter.

The measurements of DHEA-S are widely used in clinical practice. Elevated concentrations of this steroid are found in patients with adrenal hyperplasia, adrenocortical carcinoma, or hirsutism. Low levels of DHEA-S are detectable in patients suffering from adrenal dysfunction or hypopituitarism.

3. TEST PRINCIPLE

Solid phase enzyme-linked immunosorbent assay (ELISA) based on the competition principle. An unknown amount of antigen present in the sample and a fixed amount of enzyme labelled antigen compete for the binding sites of the antibodies coated onto the wells. After incubation the wells are washed to stop the competition reaction. After the substrate reaction the intensity of the developed color is inversely proportional to the amount of the antigen in the sample. Results of samples can be determined directly using the standard curve.

4. WARNINGS AND PRECAUTIONS

1. For *in-vitro diagnostic* use only. For professional use only.
2. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.
3. In case of severe damage of the kit package please contact IBL or your supplier in written form, latest one week after receiving the kit. Do not use damaged components in test runs, but keep safe for complaint related issues.
4. Obey lot number and expiry date. Do not mix reagents of different lots. Do not use expired reagents.
5. Follow good laboratory practice and safety guidelines. Wear lab coats, disposable latex gloves and protective glasses where necessary.
6. Reagents of this kit containing hazardous material may cause eye and skin irritations. See MATERIALS SUPPLIED and labels for details. Material Safety Data Sheets for this product are available on the IBL-Homepage or upon request directly from IBL.
7. Chemicals and prepared or used reagents have to be treated as hazardous waste according to national biohazard and safety guidelines or regulations.
8. The cleaning staff should be guided by the professionals regarding potential hazards and handling.
9. Avoid contact with Stop solution. It may cause skin irritations and burns.
10. All reagents of this kit containing human serum or plasma have been tested and were found negative for anti-HIV I/II, HBsAg and anti-HCV. However, a presence of these or other infectious agents cannot be excluded absolutely. For this reason reagents should be treated as potential biohazards in use and for disposal.

5. STORAGE AND STABILITY

The kit is shipped at ambient temperature and should be stored at 2-8 °C. Keep away from heat or direct sunlight. The storage and stability of specimens and prepared reagents is stated in the corresponding chapters.

The microtiter strips are stable up to the expiry date of the kit in the broken, but tightly closed bag when stored at 2-8 °C.

6. SPECIMEN COLLECTION AND STORAGE

Serum, Plasma (EDTA, Heparin, Citrate)

The usual precautions for venipuncture should be observed. It is important to preserve the chemical integrity of a blood specimen from the moment it is collected until it is assayed. Do not use grossly hemolytic, icteric or grossly lipemic specimens. Samples appearing turbid should be centrifuged before testing to remove any particulate material.

Storage:	2-8°C	≤ -20°C (Aliquots)	Keep away from heat or direct sunlight. Avoid repeated freeze-thaw cycles.
Stability:	7 days	6 months	

7. MATERIALS SUPPLIED

Quantity	Symbol	Component
1 x 12 x 8	MTP	Microtiter Plate Break apart strips. Coated with anti-DHEA-S antibody (polyclonal).
1 x 6 x 1 mL	CAL A-F	Standard A-F 0; 0.1; 0.3; 0.9; 2.7; 10.0 µg/mL Ready to use. Contains: DHEA-S, Human serum, stabilizers.
1 x 2 x 1 mL	CONTROL 1+2	Control 1+2 Ready to use. Contains: DHEA-S, low and high, Human serum, stabilizers. Exact concentrations see vial labels or QC certificate.
1 x 25 mL	ENZCONJ	Enzyme Conjugate Ready to use. Contains: DHEA-S conjugated to HRP, stabilizers.
1 x 50 mL	WASHBUF CONC	Wash Buffer Concentrate (20x) Contains: Tween, stabilizers.
1 x 15 mL	TMB SUBS	TMB Substrate Solution Ready to use. Contains: TMB, Buffer, stabilizers.
1 x 15 mL	TMB STOP	TMB Stop Solution Ready to use. 1 M H ₂ SO ₄ .
3 x	FOIL	Adhesive Foil

8. MATERIALS REQUIRED BUT NOT SUPPLIED


1. Micropipettes (Multipette Eppendorf or similar devices, < 3 % CV). Volume: 25; 100; 200 µL
2. Vortex mixer
3. 8-Channel Micropipettor with reagent reservoirs
4. Wash bottle, automated or semi-automated microtiter plate washing system
5. Microtiter plate reader capable of reading absorbance at 450 nm (reference wavelength 600-650 nm)
6. Bidistilled or deionised water
7. Paper towels, pipette tips and timer

9. PROCEDURE NOTES

1. Any improper handling of samples or modification of the test procedure may influence the results. The indicated pipetting volumes, incubation times, temperatures and pretreatment steps have to be performed strictly according to the instructions. Use calibrated pipettes and devices only.
2. Once the test has been started, all steps should be completed without interruption. Make sure that required reagents, materials and devices are prepared ready at the appropriate time. Allow all reagents and specimens to reach room temperature (18-25 °C) and gently swirl each vial of liquid reagent and sample before use. Mix reagents without foaming.
3. Avoid contamination of reagents, pipettes and wells/tubes. Use new disposable plastic pipette tips for each component and specimen. Do not interchange caps. Always cap not used vials. Do not reuse wells/tubes or reagents.
4. It is advised to determine samples in duplicate to be able to identify potential pipetting errors.
5. Use a pipetting scheme to verify an appropriate plate layout.
6. Incubation time affects results. All wells should be handled in the same order and time sequences. It is recommended to use an 8-channel Micropipettor for pipetting of solutions in all wells.

7. Microtiter plate washing is important. Improperly washed wells will give erroneous results. It is recommended to use a multichannel pipette or an automatic microtiter plate washing system. Do not allow the wells to dry between incubations. Do not scratch coated wells during rinsing and aspiration. Rinse and fill all reagents with care. While rinsing, check that all wells are filled precisely with Wash Buffer, and that there are no residues in the wells.
8. Humidity affects the coated wells/tubes. Do not open the pouch until it reaches room temperature. Unused wells/tubes should be returned immediately to the resealed pouch including the desiccant.

10. PRE-TEST SETUP INSTRUCTIONS

	The contents of the kit for 96 determinations can be divided into 3 separate runs. The volumes stated below are for one run with all strips (96 determinations).
---	--

10.1. Preparation of lyophilized or concentrated components

Dilute / dissolve	Component		Diluent	Relation	Remarks	Storage	Stability
50 mL	WASHBUF CONC	ad 1000 mL	bidist. water	1:20	Mix vigorously.	2-8°C	4 weeks

10.2. Dilution of Samples

Samples suspected to contain concentrations higher than the highest standard have to be diluted with Standard A. Dilution has to be made in glass tubes. Measured results must be multiplied with the dilution factor to obtain corrected results.

11. TEST PROCEDURE

1.	Pipette 25 µL of each Standard, Control and sample into the respective wells of the Microtiter Plate.
2.	Pipette 200 µL of Enzyme Conjugate into each well.
3.	Thoroughly mix for 10 seconds.
4.	Cover plate with adhesive foil. Incubate 60 min at RT (18-25°C).
5.	Remove adhesive foil. Discard incubation solution. Wash plate 4 x with 250 µL of diluted Wash Buffer . Remove excess solution by tapping the inverted plate on a paper towel.
6.	Pipette 100 µL of TMB Substrate Solution into each well.
7.	Incubate 15 min at RT (18-25°C) (manual). In case of automation of the assay the incubation time can be reduced to 12 min.
8.	Stop the substrate reaction by adding 100 µL of TMB Stop Solution into each well. Briefly mix contents by gently shaking the plate. Colour changes from blue to yellow.
9.	Measure optical density with a photometer at 450 nm (Reference-wavelength: 600-650 nm) within 15 min after pipetting of the Stop Solution.

12. QUALITY CONTROL

The test results are only valid if the test has been performed following the instructions. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or comparable standards/laws. User and/or laboratory must have a validated system to get diagnosis according to GLP. All kit controls must be found within the acceptable ranges as stated on the labels and the QC certificate. If the criteria are not met, the run is not valid and should be repeated. Each laboratory should use known samples as further controls. It is recommended to participate at appropriate quality assessment trials.

In case of any deviation the following technical issues should be proven: Expiration dates of (prepared) reagents, storage conditions, pipettes, devices, incubation conditions and washing methods.

13. CALCULATION OF RESULTS

The obtained OD of the standards (y-axis, linear) are plotted against their concentration (x-axis, logarithmic) either on semi-logarithmic graph paper or using an automated method. A good fit is provided with cubic spline, 4 Parameter Logistics or Logit-Log.

For the calculation of the standard curve, apply each signal of the standards (one obvious outlier of duplicates might be omitted and the more plausible single value might be used).

The concentration of the samples can be read directly from the standard curve.

Samples showing concentrations above the highest standard have to be diluted as described in PRE-TEST SETUP INSTRUCTIONS and reassayed.

In case of diluted samples the values have to be multiplied with the corresponding dilution factor.

Conversion:

DHEA-S ($\mu\text{g/mL}$) \times 2.60 = $\mu\text{mol/L}$


Reportable Ranges:

Serum: 0.05 – 10.0 $\mu\text{g/mL}$ DHEA-S

Typical Calibration Curve

(Example. Do not use for calculation!)

Standard	DHEA-S ($\mu\text{g/mL}$)	OD _{Mean}	OD/OD _{max} (%)
A	0.0	2.193	100
B	0.1	1.664	76
C	0.3	1.265	58
D	0.9	0.804	37
E	2.7	0.424	19
F	10.0	0.172	8



14. EXPECTED VALUES

The results themselves should not be the only reason for any therapeutical consequences. They have to be correlated to other clinical observations and diagnostic tests.

Apparently healthy subjects show the following values:

	n	DHEA-S [$\mu\text{g/mL}$]	
		Range	Median
Women	52	0.37 – 2.71	1.15
Men	45	0.73 – 3.81	1.90

It is recommended that each laboratory establishes its own range of normal values.

15. LIMITATIONS OF THE PROCEDURE

Specimen collection and storage have a significant effect on the test results. See SPECIMEN COLLECTION AND STORAGE for details.

For cross-reactivities, see PERFORMANCE.

Azide and Thimerosal in concentrations $>$ 0.1% interfere with the assay and cause falsely results.

The following blood components do not have a significant effect (\pm 20% of expected) on the test results up to the below stated concentrations:

Hemoglobin	2 mg/mL
Bilirubin	1 mg/mL
Triglyceride	30 mg/mL

16. PERFORMANCE


Analytical Specificity (Cross Reactivity)	Substance		Cross Reactivity (%)	Cross-reactivity of other substances tested \leq 0.01 %
	Dehydroepiandrosterone sulphate		100.0	
	Androsterone sulfate		5.67	
	17- α -Hydroxyprogesteronsulfat		0.13	
	Estrone		2.62	
	Testosterone		2.13	
Progesterone		0.93		
Analytical Sensitivity (Limit of Detection)	Serum	0.004 $\mu\text{g/mL}$	Mean signal (Zero-Standard) - 2SD	

Precision	Serum (µg/mL)			
	Mean	SD	CV (%)	N
Intra-Assay	0.89	0.06	6.6	18
	3.14	0.14	4.3	22
	9.34	0.51	5.4	18
Inter-Assay	0.9	0.05	4.9	20
	1.9	0.10	5.2	20
	3.1	0.14	4.6	20
Linearity	Serum			
	Dilution	Measured (µg/mL)	Recovery (%)	
	-	4.9	100	
	1:2	2.2	90	
	1:4	1.0	82	
	1:8	0.5	82	
	1:16	0.3	98	
	-	3.5	100	
	1:2	1.7	97	
	1:4	0.8	91	
	1:8	0.4	91	
	1:16	0.2	91	
	-	3.7	100	
	1:2	1.8	97	
	1:4	1.0	108	
1:8	0.5	108		
1:16	0.2	86		
Recovery	Conc. (µg/mL)	Added (µg/mL)	Measured (µg/mL)	Recovery (%)
	Serum 1 (0)	0.1	0.2	95
		0.64	0.6	94
		1.93	1.8	93
		7.14	6.6	92
	Serum 2 (0.9)	0.1	1.2	108
		0.64	1.6	104
		1.93	3.3	117
		7.14	8.7	108
	Serum 3 (3.7)	0.1	3.8	97
		0.64	4.5	104
		1.93	5.6	99
7.14		10.5	97	
Method Comparison	Serum	IBL-ELISA = 0.96 x DRG-ELISA + 0.21		r = 0.95; n = 77


17. PRODUCT LITERATURE REFERENCES

- Buettner R., Bollheimer L. C., Zietz B., Drobnik W., Lackner K., Schmitz G., Schölmerich J., Palitzsch K.D.: Definition and characterization of relative hypo- and hyperleptinemia in a large Caucasian population, *Journal of Endocrinology* 2002; 175(3):745–756
- Martin A., Saathoff M., Kuhn F. et al.: A Functional ABCC11 Allele Is Essential in the Biochemical Formation of Human Axillary Odor, *Journal of Investigative Dermatology* 2010, Volume 130:529–540
- Sipilä S.; Heikkinen E., Cheng S., Suominen H., Saari P., Kovanen V., Alén M., Ratanen T. Endogenous hormones, muscle strength, and risk of fall-related fractures in older women, *Journal of Gerontology: Medical Sciences* 2006, Vol. 61A, No. 1:92–96
- Thomas Lothar; *Labor und Diagnose: Indikation und Bewertung von Laborbefunden für die medizinische Diagnostik*, 7. Auflage, 2008
- Tietz, N.W. *Textbook of Clinical Chemistry*. Saunders, 1986.
- Weidler C., Struharova S., Schmidt M., Ugele B., Schölmerich J., Straub R. H.: Tumor Necrosis Factor Inhibits Conversion of Dehydroepiandrosterone Sulfate (DHEAS) to DHEA in Rheumatoid Arthritis Synovial Cells. A Prerequisite for Local Androgen Deficiency, *Arthritis & Rheumatism* 2005, Vol. 52, No. 6:1721–1729

Symbols / Symbole / Symbôles / Símbolos / Símbolos / Σύμβολα

	Cat.-No.: / Kat.-Nr.: / No.- Cat.: / Cat.-No.: / N.º Cat.: / N.-Cat.: / Αριθμός-Κατ.:
	Lot-No.: / Chargen-Bez.: / No. Lot: / Lot-No.: / Lote N.º: / Lotto n.: / Αριθμός -Παραγωγή:
	Use by: / Verwendbar bis: / Utiliser à: / Usado por: / Usar até: / Da utilizzare entro: / Χρησιμοποιείται από:
	No. of Tests: / Kitgröße: / Nb. de Tests: / No. de Determ.: / N.º de Testes: / Quantità dei tests: / Αριθμός εξετάσεων:
	Concentrate / Konzentrat / Concentré / Concentrar / Concentrado / Concentrato / Συμπύκνωμα
	Lyophilized / Lyophilisat / Lyophilisé / Liofilizado / Liofilizado / Liofilizzato / Λυοφιλιασμένο
	In Vitro Diagnostic Medical Device. / In-vitro-Diagnostikum. / Appareil Médical pour Diagnostics In Vitro. / Dispositivo Médico para Diagnóstico In Vitro. / Equipamento Médico de Diagnóstico In Vitro. / Dispositivo Medico Diagnostico In vitro. / Ιατρική συσκευή για In-Vitro Διάγνωση.
	Evaluation kit. / Nur für Leistungsbewertungszwecke. / Kit pour évaluation. / Juego de Reactivos para Evaluació. / Kit de avaliação. / Kit di valutazione. / Κιτ Αξιολόγησης.
	Read instructions before use. / Arbeitsanleitung lesen. / Lire la fiche technique avant emploi. / Lea las instrucciones antes de usar. / Ler as instruções antes de usar. / Leggere le istruzioni prima dell'uso. / Διαβάστε τις οδηγίες πριν την χρήση.
	Keep away from heat or direct sun light. / Vor Hitze und direkter Sonneneinstrahlung schützen. / Garder à l'abri de la chaleur et de toute exposition lumineuse. / Manténgase alejado del calor o la luz solar directa. / Manter longe do calor ou luz solar directa. / Non esporre ai raggi solari. / Να φυλάσσεται μακριά από θερμότητα και άμεση επαφή με το φως του ηλίου.
	Store at: / Lagern bei: / Stocker à: / Almacene a: / Armazenar a: / Conservare a: / Αποθήκευση στους:
	Manufacturer: / Hersteller: / Fabricant: / Productor: / Fabricante: / Fabbricante: / Παραγωγός:
	Caution! / Vorsicht! / Attention! / ¡Precaución! / Cuidado! / Attenzione! / Προσοχή!
<p>Symbols of the kit components see MATERIALS SUPPLIED. Die Symbole der Komponenten sind im Kapitel KOMPONENTEN DES KITS beschrieben. Voir MATERIEL FOURNI pour les symbôles des composants du kit. Símbolos de los componentes del juego de reactivos, vea MATERIALES SUMINISTRADOS. Para símbolos dos componentes do kit ver MATERIAIS FORNECIDOS. Per i simboli dei componenti del kit si veda COMPONENTI DEL KIT. Για τα σύμβολα των συστατικών του κιτ συμβουλευτείτε το ΠΑΡΕΧΟΜΕΝΑ ΥΛΙΚΑ.</p>	

IBL AFFILIATES WORLDWIDE

	IBL International GmbH Flughafenstr. 52A, 22335 Hamburg, Germany	Tel.: + 49 (0) 40 532891 -0 Fax: -11 E-MAIL: IBL@IBL-International.com WEB: http://www.IBL-International.com
	IBL International Corp. 194 Wildcat Road, Toronto, Ontario M3J 2N5, Canada	Tel.: +1 (416) 645 -1703 Fax: -1704 E-MAIL: Sales@IBL-International.com WEB: http://www.IBL-International.com

LIABILITY: Complaints will be accepted in each mode –written or vocal. Preferred is that the complaint is accompanied with the test performance and results. Any modification of the test procedure or exchange or mixing of components of different lots could negatively affect the results. These cases invalidate any claim for replacement. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the kit during transportation is not subject to the liability of the manufacturer