

— BiOLiS 50i — Superior

Automated Clinical Analyzer



TOKYO BOEKI MEDISYS

LABORATORY SOLUTION

BiOLIS 50i Superior

*Higher test performance and operability with new functions
The clinical analyzer in the next stage*





The best to the requirements for functionality and usability

BiOLiS 50i Superior

The highest throughput in BIOLIS series

- 480 tests/hour for photometry. 580 tests/hour with ISE.

Rapid correspondence by a new function for HbA1c

- Hemolysis of whole blood sample for HbA1c is available automatically.
※The throughput of HbA1c is different from that of general chemistries.

Automatic sample clot detection & automatic cleaning

- Clot detector for the sample probe enhances the accuracy of the test results.
- Automatically cleaning the sample probe in instances where clots are detected.
- Sampling ceases automatically in case clots are not removed.

Accommodate multiple languages

- Multilingual support software provides a user-friendly interface and higher usability.

User interface

- Easy touch screen operation.
- The interface includes improved and additional functions such as sample status screen, and others.

Dedicated sample probe for ISE

- The sample probe **directly connected** to ISE transfers the sample to electrodes.

Common reagent bottles

- R1: 70 ml, R2: 20 ml, two sizes can be used.


Accuracy Features

- Sample clot detection
- Dedicated sample probe for ISE
- Carry-over protection program for reaction cuvettes and probes
- Pre-dilution
- Auto-rerun
- Cuvette skip
- Five mixing speeds

Operation Features

- Sample cups (0.5 ml, 1.5 ml),
Primary tubes (2 ml, 5 ml, 7 ml, 10 ml)
- Availability for using multiple bottles of the same item
- Reaction waste stored in a dedicated tank
- Common reagent bottles (R1: 70 ml, R2: 20 ml)
- Hemolysis of whole blood sample for HbA1c
- Sample and reagent barcode readers as standard feature
- Bidirectional communication

Main Unit Arrangement



1 Reagent tray

R1: 36, R2: 35 position.
Cooling temperature: 8-12°C

2 R1 probe

20-240 μl (in 0.5 μl steps) dispensing
※Minimum total reaction volume is 120 μl

3 R2 probe

4 Reaction tray

90 semi-disposable cuvettes.
Cuvette skip function is available.

5 Cuvette washing station

Cuvette washing with alkaline and acid washing solutions and heated water.

6 Mixing unit

Five mixing speeds

7 Sample probe

1-20 μl (in 0.1 μl steps) dispensing

8 ISE module

Na, K, Cl measurement

9 ISE probe

Dedicated sample probe for ISE

10 Sample tray

72 samples on board.

Specifications

Analysis	System	Discrete single line random access multi-test analysis
	Number of test items on board	36 + 3 items (ISE)
	Throughput	480 tests/hour, 580 tests/hour with ISE *160 tests/hour for HbA1c only (with hemolysis)
	Analysis methods	End point assay, Rate assay, ISE
	Calibration curve	8 kinds (Linear, Spline, etc)
Sample	Contents of test	Serum, Plasma, Urine, Dialysis, CSF (ISE not available for CSF)
	Sample container	Sample cups (0.5 ml, 1.5 ml), Primary tubes (2ml, 5 ml, 7 ml, 10 ml)
	Number of samples on board	72 patient samples and standard samples, 30 positions for QC and/or STAT samples
	Sample dispensing volume	1- 20 µl (in 0.1 µl steps)
	Dilution ratio	3-100 times
	STAT	Available during measurement
	Sample barcode reader	Available
Reagent	Number of bottles on board	R1: 36 positions, R2: 35 positions
	Bottle size	R1: 70 ml, R2: 20 ml
	Reagent dispensing volume	20-240 µl (in 0.5 µl steps)
	Reagent Storage	Reagent tray cooling available
	Reagent residual volume	Level sensing or count down calculation
	Reagent barcode reader	Available
Reaction	Cuvette material	Plastics (Semi-disposable)
	Number of cuvettes	90
	Cuvette washing	Auto washing with heated water and 2 kinds of washing solutions
	Washing solution	Alkaline and acid washing solutions
	Cuvette skip function	Available
	Reaction time	3.5 min. after sample dispensing + 5 min. after R2 dispensing
	Reaction volume	120- 300 µl
	Reaction temperature	37°C ± 0.3°C
	Photometric assay method	Grating method
	Optical measurements	12 fixed wavelengths (340, 380, 405, 450, 505, 546, 570, 600, 660, 700, 750, 800 nm)
	Optical source	Tungsten halogen lamp
	Optical range	OD 0- 3.0
	Mixing	By stirrers
	Mixing speed	5 speed levels available
Interface	Operation	Personal computer
	OS	Windows 7
	Output	RS232C serial cable, LAN (TCP/IP)
	Reaction curve monitor	Optical absorbance graphic display
	Quality control	Realtime QC, Within a day, Day-to-day variation
	Test results storage	100,000 tests

*The specifications are subject to change without notice.

Test Items List

Clinical Chemistry	LD(LDH)	AST(GOT)	ALT(GPT)	ALP	γ-GTP	CK(CPK)	LAP
	ChE	AMY	P-AMY	BUN	CRE	UA	T-CHO
	TG	HDL-C	LDL-C	TP	ALB	Ca	IP
	Mg	Fe	UIBC	D-BIL	T-BIL	GLU	NEFA
	PL	ALD	SIA	TBA	GA	CK-MB	
	Apo A-I	Apo A-II	Apo B	Apo CII	Apo CIII	Apo E	
	Lp(a)	Fer	µTP	NAG	µALB	HbA1c	HBAO
Immunoassay	CRP	RF	ASO	TPLA	RPR	β2-m	Tf
	C3	C4	Ig-G	Ig-A	Ig-M	Ig-E	
TDM	BRP	CBZ	CSA	DIG	DIGT	HAL	PB
	PHT	THEO	VPA				
Coagulation	ATIII	APL	PLG	PC	FDP	D dimer	
ISE	Na	K	Cl				

*Some test items are during verification.

User Interface


Simple Operation

Direct access to each function by clicking tabs and buttons shown on the main screen


New feature for time management

The list of time records for orders or test results is displayed in "Sample Status Screen"

Order Entry Screen



Sample Status Screen



Test Result Screen




Run Monitor Screen



QC Graph Screen



HbA1c sample preparation and measurement




Installation

Conditions

Item	Description
Dimensions and weight	Analyzer W1050 x D750 x H1145 (mm) Approx. 300kg
Power supply	AC100/115/230V±10%, 50/60Hz Voltage fluctuation: Less than 10%
Power consumption	1.2kVA
Grounding	Earth resistance of ground terminal should be less than 10Ω.
Ambient temperature	15-30°C ±2°C/hour during measurement
Humidity	45-85% (No condensation)
Water consumption	Max. 13 ℓ/hour
Waste liquids	Separate drainage (low and high density waste)


The space required for installation



※At least 500 mm space required for safety and maintenance.

ISE module (Option)

- Direct method
- Dedicated sample probe for ISE
- Na, K, Cl in serum, plasma and urine can be measured
- Sample volume: 100 μℓ
- Electrodes can be commonly used in BIOLIS series





GLOBAL PARTNER FOR LABORATORY SOLUTION

ISO9001: 2008
ISO13485: 2003



REPRESENTED BY
TOKYO BOEKI MACHINERY LTD.
MEDICAL INSTRUMENTS DIVISION
Headquarters
2-13-8, Hatchobori, Chuo-ku, Tokyo 104-8510, Japan
Tel: +81-3-3555-7353 Fax: +81-3-3555-7367
URL <http://www.tokyo-boeki-machinery.jp>



DEVELOPED AND MANUFACTURED BY
TOKYO BOEKI MEDISYS INC.
Headquarters
1-14-21, Higashi-Toyoda, Hino, Tokyo 191-0052, Japan
Tel: +81-42-587-7777 Fax: +81-42-587-7781
URL http://www.tb-medisys.co.jp/index_en.html

