

## Instrument information

Features	Available accessories
Colour touch display	QuikRead go Accumulator
Result memory for 100 patient results and 100 QC results	Barcode reader
POCT1-A2 bidirectional connectivity LIS01-A2 unidirectional connectivity	Printer
Size 27*15.5*14.5 cm	Workstation for QuikRead CRP
Weight 1.7 kg	QuikRead go software USB stick
	QuikRead go WLAN adapter

## Test kit information

Test	Sample type	Measuring range
QuikRead go CRP	Whole blood, serum, plasma	5–200 mg/l
QuikRead go CRP+Hb	Whole blood Serum and plasma (only CRP)	5–200 mg/l CRP 50–245 g/l Hb
QuikRead go Strep A	Throat swab	7x10 <sup>4</sup> CFU*/swab
QuikRead go iFOBT	Human faeces	Positive above 75 ng/ml (15 µg/g faeces) of which quantitative range: 75–1000 ng/ml (15–200 µg/g faeces)

\**Streptococcus pyogenes*

Orion Diagnostica develops, manufactures and markets diagnostic test systems for healthcare professionals with over 40 years of experience. Building well-being is the core of our mission. QuikRead go is a registered trademark of Orion Diagnostica Oy.



Orion Diagnostica Oy  
Koivu-Mankkaan tie 6 B, P.O. Box 83, FI-02101 Espoo, Finland  
Tel. +358 10 4261, Fax +358 10 426 2794  
www.oriondiagnostica.com, www.quikread.com

8015-06GB, 02/2016. International version. Not intended for US market.

QuikRead go  
Products



Your support in treatment decisions

# QuikRead go<sup>®</sup> Products



# Your support in treatment decisions

## *Reliable with high quality*

High quality QuikRead go test kits with a fully automatic testing procedure provide reliable, accurate results. You can trust that the results are comparable to results provided by laboratory systems and reference methods. The user interface offers easy patient and operator identification, providing further assurance. Continuous development keeps the instrument up to date for your everyday needs.

## *Fast and easy to use*

QuikRead go provides you with results within a few minutes. The fast results enable immediate diagnosis and treatment decisions during patient consultation. The interactive user interface, available in 22 languages, guides you through the easy measurement procedure. No special laboratory skills are required and the instrument does not need any regular maintenance. The QuikRead go system provides flexibility for your everyday work. When near patient results are required, you can easily carry the instrument to the point-of-care.

## *Connectable*

The system offers you a selection of connectivity features. QuikRead go instrument is connectable to most hospital (HIS) and laboratory information (LIS) systems. Bidirectional communication offers new possibilities for user, instrument and quality management. Results can be printed on labels or stored on USB sticks for future use as well.

## *Test selection*

The same QuikRead go instrument can be used with all the available QuikRead go tests and the instrument recognises automatically the test in use. The selection currently contains the following tests:

- » QuikRead go CRP
- » QuikRead go CRP+Hb
- » QuikRead go Strep A
- » QuikRead go iFOBT



QuikRead go, the easy to use point-of-care test system, brings reliability and speed to your everyday working processes. This flexible system supports you in making diagnoses and treatment decisions efficiently thus increasing patient comfort and satisfaction.